

## Process Report

**Order ID :** Q3320

| Date       | Bottle ID | Lab Sample Number | InComing Vacuum | Canister ID | Can Size | Reg ID | Reg Size | Batch Code | Certification ID                                                                                                     |
|------------|-----------|-------------------|-----------------|-------------|----------|--------|----------|------------|----------------------------------------------------------------------------------------------------------------------|
| 10/10/2025 | B2509079  | Q3320-02          | -5.9            | 10332       | 6 L      | 10767  |          | C2508001   | <b>VL042828.D</b><br>Seq :VL080725<br>TUNE : VL042823.D<br>ICAL : VL042752.D<br>CCAL : VL042824.D<br>MB : VL042825.D |
| 10/10/2025 | B2509079  | Q3320-12          | -5.1            | 10333       | 6 L      | 10528  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-13          | -4.3            | 10609       | 6 L      | 10707  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-11          | -5.1            | 10607       | 6 L      | 10617  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-10          | -5.1            | 10604       | 6 L      | 10648  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-03          | -5.1            | 10602       | 6 L      | 10109  |          | C2509002   | <b>VL042942.D</b><br>Seq :VL091825<br>TUNE : VL042938.D<br>ICAL : VL042934.D<br>CCAL : VL042939.D<br>MB : VL042940.D |

## Process Report

**Order ID :** Q3320

| Date       | Bottle ID | Lab Sample Number | InComing Vacuum | Canister ID | Can Size | Reg ID | Reg Size | Batch Code | Certification ID                                                                                                     |
|------------|-----------|-------------------|-----------------|-------------|----------|--------|----------|------------|----------------------------------------------------------------------------------------------------------------------|
| 10/10/2025 | B2509079  | Q3320-09          | -5.1            | 10601       | 6 L      | 10550  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-08          | -5.5            | 10588       | 6 L      | 10525  |          | C2506003   | <b>VL042702.D</b><br>Seq :VL070325<br>TUNE : VL042698.D<br>ICAL : VL042694.D<br>CCAL : VL042699.D<br>MB : VL042700.D |
| 10/10/2025 | B2509079  | Q3320-04          | -5.1            | 10445       | 6 L      | 10185  |          | C2508001   | <b>VL042828.D</b><br>Seq :VL080725<br>TUNE : VL042823.D<br>ICAL : VL042752.D<br>CCAL : VL042824.D<br>MB : VL042825.D |
| 10/10/2025 | B2509079  | Q3320-06          | -3.9            | 10443       | 6 L      | 10613  |          | C2509002   | <b>VL042942.D</b><br>Seq :VL091825<br>TUNE : VL042938.D<br>ICAL : VL042934.D<br>CCAL : VL042939.D<br>MB : VL042940.D |
| 10/10/2025 | B2509079  | Q3320-07          | -4.3            | 10311       | 6 L      | 10226  |          | C2507001   | <b>VL042703.D</b><br>Seq :VL070325<br>TUNE : VL042698.D<br>ICAL : VL042694.D<br>CCAL : VL042699.D<br>MB : VL042700.D |
| 10/10/2025 | B2509079  | Q3320-05          | -5.1            | 10285       | 6 L      | 10579  |          | C2509001   | <b>VL042907.D</b><br>Seq :VL090825<br>TUNE : VL042899.D<br>ICAL : VL042847.D<br>CCAL : VL042900.D<br>MB : VL042901.D |
| 10/10/2025 | B2509079  | Q3320-01          | -4.3            | 10053       | 6 L      | 10616  |          | C2509002   | <b>VL042942.D</b><br>Seq :VL091825<br>TUNE : VL042938.D<br>ICAL : VL042934.D<br>CCAL : VL042939.D<br>MB : VL042940.D |