

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL100225

Review By	Semsettin Yesilyurt	Review On	10/3/2025 9:49:03 AM
Supervise By	Maresh Dadoda	Supervise On	10/10/2025 2:54:57 AM
SubDirectory	VL100225	HP Acquire Method	TO15AIR
		HP Processing Method	VL100225AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2682
Initial Calibration Stds	AP2678,AP2679,AP2680
CCC	AP2678
Internal Standard/PEM	AP2682
ICV/I.BLK	AP2681
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	Sampled	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL043030.D	02 Oct 2025 09:11		SY/MD	Ok
2	VSTDICC0.5	VSTDICC0.5	VL043031.D	02 Oct 2025 09:50		SY/MD	Ok,M
3	VSTDICCC010	VSTDICCC010	VL043032.D	02 Oct 2025 10:23	LR-13,80 and QR-79	SY/MD	Ok,M
4	VSTDICC002	VSTDICC002	VL043033.D	02 Oct 2025 10:57		SY/MD	Ok,M
5	VSTDICC001	VSTDICC001	VL043034.D	02 Oct 2025 11:29		SY/MD	Ok,M
6	VSTDICC0.5	VSTDICC0.5	VL043035.D	02 Oct 2025 12:02	Not used	SY/MD	Not Ok
7	VSTDICC0.1	VSTDICC0.1	VL043036.D	02 Oct 2025 12:35		SY/MD	Ok,M
8	VSTDICC0.03	VSTDICC0.03	VL043037.D	02 Oct 2025 13:07		SY/MD	Ok,M
9	VSTDICC015	VSTDICC015	VL043038.D	02 Oct 2025 13:41		SY/MD	Ok,M
10	VSTDICV010	ICVVL100225	VL043039.D	02 Oct 2025 14:14		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL100825

Review By	Semsettin Yesilyurt	Review On	10/9/2025 1:56:13 PM		
Supervise By	Mahesh Dadoda	Supervise On	10/10/2025 4:42:55 AM		
SubDirectory	VL100825	HP Acquire Method	TO15AIR	HP Processing Method	VL100225AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2682
Initial Calibration Stds	AP2678,AP2679,AP2680
CCC	AP2678
Internal Standard/PEM	AP2682
ICV/I.BLK	AP2681
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL043053.D	08 Oct 2025 08:00		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL043054.D	08 Oct 2025 08:54		SY/MD	Ok,M
3	VL1008ABL01	VL1008ABL01	VL043055.D	08 Oct 2025 09:39		SY/MD	Ok
4	VL1008ABS01	VL1008ABS01	VL043056.D	08 Oct 2025 10:23		SY/MD	Ok,M
5	QC100625 CAN 10595	QC100625 CAN 10595	VL043057.D	08 Oct 2025 11:10	lab contaminate	SY/MD	Not Ok
6	QC100725 CAN 10420	QC100725 CAN 10420	VL043058.D	08 Oct 2025 11:43		SY/MD	Ok
7	QC100725 CAN 10133	QC100725 CAN 10133	VL043059.D	08 Oct 2025 12:16		SY/MD	Ok
8	QC100725 CAN 10122	QC100725 CAN 10122	VL043060.D	08 Oct 2025 12:49		SY/MD	Ok
9	Q3306-01	IA1	VL043061.D	08 Oct 2025 13:23		SY/MD	Ok,M
10	Q3306-01DUP	IA1DUP	VL043062.D	08 Oct 2025 13:58		SY/MD	Ok,M
11	Q3306-02	OA1	VL043063.D	08 Oct 2025 14:32		SY/MD	Ok,M
12	Q3306-04	SV2	VL043064.D	08 Oct 2025 15:05	rerun lower dilution	SY/MD	Not Ok
13	Q3306-03	SV1	VL043065.D	08 Oct 2025 15:38	rerun lower dilution	SY/MD	Not Ok
14	Q3306-04	SV2	VL043066.D	08 Oct 2025 19:57		SY/MD	Ok,M
15	Q3306-03	SV1	VL043067.D	08 Oct 2025 20:30		SY/MD	Ok,M
16	VIBLK	VIBLK	VL043068.D	08 Oct 2025 21:03		SY/MD	Ok
17	QC100684 CAN 10595	QC100684 CAN 10595	VL043069.D	08 Oct 2025 21:36	lab contaminate	SY/MD	Not Ok
18	Q3320-13	OA1	VL043070.D	08 Oct 2025 22:10		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL100825

Review By	Semsettin Yesilyurt	Review On	10/9/2025 1:56:13 PM		
Supervise By	Mahesh Dadoda	Supervise On	10/10/2025 4:42:55 AM		
SubDirectory	VL100825	HP Acquire Method	TO15AIR	HP Processing Method	VL100225AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2682				
Initial Calibration Stds	AP2678,AP2679,AP2680				
CCC	AP2678				
Internal Standard/PEM	AP2682				
ICV/I.BLK	AP2681				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q3320-13DL	OA1DL	VL043071.D	08 Oct 2025 22:43	Not required	SY/MD	Not Ok
20	Q3320-09	IA5	VL043072.D	08 Oct 2025 23:16		SY/MD	Ok,M
21	Q3320-09DL	IA5DL	VL043073.D	08 Oct 2025 23:49	Not required	SY/MD	Not Ok
22	Q3320-05	IA3	VL043074.D	09 Oct 2025 00:22		SY/MD	Ok,M
23	Q3320-05DL	IA3DL	VL043075.D	09 Oct 2025 00:55	Not required	SY/MD	Not Ok
24	Q3320-07	IA4	VL043076.D	09 Oct 2025 01:28		SY/MD	Ok,M
25	Q3320-07DL	IA4DL	VL043077.D	09 Oct 2025 02:01	Not required	SY/MD	Not Ok

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL100925

Review By	Semsettin Yesilyurt	Review On	10/10/2025 2:03:15 PM
Supervise By	Maresh Dadoda	Supervise On	10/10/2025 11:17:18 PM
SubDirectory	VL100925	HP Acquire Method	TO15AIR
		HP Processing Method	VL100225AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2683
Initial Calibration Stds	AP2678,AP2679,AP2680
CCC	AP2678
Internal Standard/PEM	AP2683
ICV/I.BLK	AP2681
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL043078.D	09 Oct 2025 07:56		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL043079.D	09 Oct 2025 09:29		SY/MD	Ok,M
3	VL1009ABL01	VL1009ABL01	VL043080.D	09 Oct 2025 10:10		SY/MD	Ok
4	VL1009ABS01	VL1009ABS01	VL043081.D	09 Oct 2025 10:54		SY/MD	Ok,M
5	QC100625 CAN 10595	QC100625 CAN 10595	VL043082.D	09 Oct 2025 11:35		SY/MD	Ok
6	Q3320-01	IA1	VL043083.D	09 Oct 2025 12:09		SY/MD	Ok,M
7	Q3320-01DUP	IA1DUP	VL043084.D	09 Oct 2025 12:44		SY/MD	Ok,M
8	Q3320-03	IA2	VL043085.D	09 Oct 2025 13:18		SY/MD	Ok,M
9	Q3320-11	IA6	VL043086.D	09 Oct 2025 13:52		SY/MD	Ok,M
10	Q3320-04	SV2	VL043087.D	09 Oct 2025 14:25	Need 40X	SY/MD	Dilution
11	Q3320-06	SV3	VL043088.D	09 Oct 2025 14:57	Run @ Lower Dilution	SY/MD	Not Ok
12	Q3320-02	SV1	VL043089.D	09 Oct 2025 15:31	Need 200X	SY/MD	Dilution
13	Q3320-08	SV4	VL043090.D	09 Oct 2025 16:07	Need further dilution	SY/MD	Not Ok
14	Q3320-10	SV5	VL043091.D	09 Oct 2025 16:42	Run @ Lower Dilution	SY/MD	Not Ok
15	Q3320-12DL	SV6DL	VL043092.D	09 Oct 2025 17:17		SY/MD	Ok,M
16	Q3320-04DL	SV2DL	VL043093.D	09 Oct 2025 18:20		SY/MD	Ok
17	Q3320-06	SV3	VL043094.D	09 Oct 2025 18:54		SY/MD	Ok,M
18	Q3320-08	SV4	VL043095.D	09 Oct 2025 19:29	Need 40X	SY/MD	Dilution

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL100925

Review By	Semsettin Yesilyurt	Review On	10/10/2025 2:03:15 PM		
Supervise By	Mahesh Dadoda	Supervise On	10/10/2025 11:17:18 PM		
SubDirectory	VL100925	HP Acquire Method	TO15AIR	HP Processing Method	VL100225AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2683				
Initial Calibration Stds	AP2678,AP2679,AP2680				
CCC	AP2678				
Internal Standard/PEM	AP2683				
ICV/I.BLK	AP2681				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q3320-10	SV5	VL043096.D	09 Oct 2025 20:04		SY/MD	Ok,M
20	Q3320-12	SV6	VL043097.D	09 Oct 2025 20:39	Need 10X	SY/MD	Dilution
21	Q3320-02DL	SV1DL	VL043098.D	09 Oct 2025 21:14	Need further dilution	SY/MD	Not Ok
22	Q3320-02DL	SV1DL	VL043099.D	09 Oct 2025 21:49	Not match	SY/MD	Not Ok
23	VIBLK	VIBLK	VL043100.D	09 Oct 2025 22:23		SY/MD	Ok
24	Q3320-08DL	SV4DL	VL043101.D	09 Oct 2025 22:58		SY/MD	Ok,M
25	Q3320-02DL	SV1DL	VL043102.D	10 Oct 2025 07:46		SY/MD	Ok

M : Manual Integration