

RETENTION TIMES OF INITIAL CALIBRATION

Lab Name:	<u>Alliance</u>	Contract:	<u>ARAM01</u>
Lab Code:	<u>ACE</u>	SDG NO.:	<u>Q3368</u>
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):	<u>11/03/2025</u> <u>11/03/2025</u>
		Calibration Times:	<u>09:40</u> <u>11:17</u>

GC Column: RTX-CLP ID: 0.32 (mm)

LAB FILE ID:	RT 200 = <u>PS032226.D</u>	RT 500 = <u>PS032227.D</u>
RT 750 = <u>PS032228.D</u>	RT 1000 = <u>PS032229.D</u>	RT 1500 = <u>PS032230.D</u>

COMPOUND	RT 200	RT 500	RT 750	RT 1000	RT 1500	MEAN RT	RT WINDOW	
							FROM	TO
2,4,5-TP(Silvex)	9.28	9.28	9.28	9.28	9.28	9.28	9.18	9.38
2,4-D	8.43	8.42	8.41	8.41	8.41	8.42	8.32	8.52
2,4-DCAA	7.29	7.29	7.28	7.28	7.28	7.28	7.18	7.38

RETENTION TIMES OF INITIAL CALIBRATION

Lab Name:	<u>Alliance</u>	Contract:	<u>ARAM01</u>
Lab Code:	<u>ACE</u>	SDG NO.:	<u>Q3368</u>
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):	<u>11/03/2025</u> <u>11/03/2025</u>
		Calibration Times:	<u>09:40</u> <u>11:17</u>

GC Column: RTX-CLP2 ID: 0.32 (mm)

LAB FILE ID:	RT 200 = <u>PS032226.D</u>	RT 500 = <u>PS032227.D</u>
RT 750 = <u>PS032228.D</u>	RT 1000 = <u>PS032229.D</u>	RT 1500 = <u>PS032230.D</u>

COMPOUND	RT 200	RT 500	RT 750	RT 1000	RT 1500	MEAN RT	RT WINDOW	
							FROM	TO
2,4,5-TP(Silvex)	9.75	9.75	9.74	9.74	9.74	9.75	9.65	9.85
2,4-D	8.89	8.87	8.86	8.86	8.86	8.87	8.77	8.97
2,4-DCAA	7.63	7.62	7.62	7.62	7.61	7.62	7.52	7.72