

Surrogate Summary

SDG No.: Q3368

Client: Aramark Uniforms

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS032225.D	PIBLK-PS032225.D	2,4-DCAA	1	500	359	72		61	136
		2,4-DCAA	2	500	508	102		61	136
I.BLK-PS032266.D	PIBLK-PS032266.D	2,4-DCAA	1	500	361	72		61	136
		2,4-DCAA	2	500	527	105		61	136
PB170386BS	PB170386BS	2,4-DCAA	1	500	556	111		61	136
		2,4-DCAA	2	500	561	112		61	136
Q3368-06	OUTSIDE-DUMPSTER	2,4-DCAA	1	500	347	69		61	136
		2,4-DCAA	2	500	439	88		61	136
Q3368-06MS	OUTSIDE-DUMPSTERMS	2,4-DCAA	1	500	400	80		61	136
		2,4-DCAA	2	500	496	99		61	136
Q3368-06MSD	OUTSIDE-DUMPSTERMSD	2,4-DCAA	1	500	399	80		61	136
		2,4-DCAA	2	500	470	94		61	136
Q3368-07	PRESS-SLUDGE	2,4-DCAA	1	500	415	83		61	136
		2,4-DCAA	2	500	478	96		61	136
I.BLK-PS032275.D	PIBLK-PS032275.D	2,4-DCAA	1	500	407	81		61	136
		2,4-DCAA	2	500	489	98		61	136
I.BLK-PS032278.D	PIBLK-PS032278.D	2,4-DCAA	1	500	412	82		61	136
		2,4-DCAA	2	500	531	106		61	136
PB170386BL	PB170386BL	2,4-DCAA	1	500	377	75		61	136
		2,4-DCAA	2	500	500	100		61	136
PB170172TB	PB170172TB	2,4-DCAA	1	500	374	75		61	136
		2,4-DCAA	2	500	498	100		61	136
I.BLK-PS032287.D	PIBLK-PS032287.D	2,4-DCAA	1	500	389	78		61	136
		2,4-DCAA	2	500	520	104		61	136