

Surrogate Summary

SDG No.: **Q3368**

Client: **Aramark Uniforms**

Analytical Method: **8081B**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PD090647.D	PIBLK-PD090647.D	Decachlorobiphen	1	20	21.3	106		57	171
		Tetrachloro-m-xyl	1	20	19.4	97		61	148
		Decachlorobiphen	2	20	20.3	102		57	171
		Tetrachloro-m-xyl	2	20	20.0	100		61	148
I.BLK-PD090935.D	PIBLK-PD090935.D	Decachlorobiphen	1	20	21.8	109		57	171
		Tetrachloro-m-xyl	1	20	19.3	96		61	148
		Decachlorobiphen	2	20	26.1	130		57	171
		Tetrachloro-m-xyl	2	20	21.2	106		61	148
Q3368-06	OUTSIDE-DUMPSTER	Decachlorobiphen	1	20	7.53	38	*	57	171
		Tetrachloro-m-xyl	1	20	9.18	46	*	61	148
		Decachlorobiphen	2	20	4.03	20	*	57	171
		Tetrachloro-m-xyl	2	20	11.0	55	*	61	148
Q3368-06MS	OUTSIDE-DUMPSTERMS	Decachlorobiphen	1	20	9.42	47	*	57	171
		Tetrachloro-m-xyl	1	20	12.8	64		61	148
		Decachlorobiphen	2	20	6.23	31	*	57	171
		Tetrachloro-m-xyl	2	20	14.7	73		61	148
Q3368-06MSD	OUTSIDE-DUMPSTERMSD	Decachlorobiphen	1	20	9.69	48	*	57	171
		Tetrachloro-m-xyl	1	20	13.3	67		61	148
		Decachlorobiphen	2	20	7.01	35	*	57	171
		Tetrachloro-m-xyl	2	20	15.1	76		61	148
Q3368-07	PRESS-SLUDGE	Decachlorobiphen	1	20	5.30	27	*	57	171
		Tetrachloro-m-xyl	1	20	3.01	15	*	61	148
		Decachlorobiphen	2	20	4.18	21	*	57	171
		Tetrachloro-m-xyl	2	20	7.46	37	*	61	148
I.BLK-PD090948.D	PIBLK-PD090948.D	Decachlorobiphen	1	20	17.7	89		57	171
		Tetrachloro-m-xyl	1	20	19.2	96		61	148
		Decachlorobiphen	2	20	12.9	65		57	171
		Tetrachloro-m-xyl	2	20	19.9	99		61	148
PB170387BS	PB170387BS	Decachlorobiphen	1	20	25.6	128		57	171
		Tetrachloro-m-xyl	1	20	26.4	132		61	148
		Decachlorobiphen	2	20	23.9	119		57	171
		Tetrachloro-m-xyl	2	20	28.0	140		61	148
PB170387BL	PB170387BL	Decachlorobiphen	1	20	18.9	95		57	171
		Tetrachloro-m-xyl	1	20	18.9	95		61	148
		Decachlorobiphen	2	20	17.1	86		57	171
		Tetrachloro-m-xyl	2	20	19.9	100		61	148
PB170172TB	PB170172TB	Decachlorobiphen	1	20	18.6	93		57	171
		Tetrachloro-m-xyl	1	20	19.3	97		61	148
		Decachlorobiphen	2	20	17.3	86		57	171
		Tetrachloro-m-xyl	2	20	19.8	99		61	148
I.BLK-PD090956.D	PIBLK-PD090956.D	Decachlorobiphen	1	20	19.2	96		57	171

Surrogate Summary

SDG No.: Q3368

Client: Aramark Uniforms

Analytical Method: 8081B

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PD090956.D	PIBLK-PD090956.D	Tetrachloro-m-xyl	1	20	19.6	98		61	148
		Decachlorobiphen	2	20	17.8	89		57	171
I.BLK-PD090959.D	PIBLK-PD090959.D	Tetrachloro-m-xyl	2	20	20.0	100		61	148
		Decachlorobiphen	1	20	19.7	98		57	171
		Tetrachloro-m-xyl	1	20	19.9	100		61	148
		Decachlorobiphen	2	20	20.1	101		57	171
		Tetrachloro-m-xyl	2	20	21.1	106		61	148
Q3368-06RE	OUTSIDE-DUMPSTERRE	Decachlorobiphen	1	20	8.73	44	*	57	171
		Tetrachloro-m-xyl	1	20	9.68	48	*	61	148
		Decachlorobiphen	2	20	8.26	41	*	57	171
		Tetrachloro-m-xyl	2	20	11.8	59	*	61	148
Q3368-07RE	PRESS-SLUDGERE	Decachlorobiphen	1	20	6.41	32	*	57	171
		Tetrachloro-m-xyl	1	20	3.03	15	*	61	148
		Decachlorobiphen	2	20	5.56	28	*	57	171
		Tetrachloro-m-xyl	2	20	7.85	39	*	61	148
I.BLK-PD090968.D	PIBLK-PD090968.D	Decachlorobiphen	1	20	19.3	96		57	171
		Tetrachloro-m-xyl	1	20	19.9	100		61	148
		Decachlorobiphen	2	20	19.0	95		57	171
		Tetrachloro-m-xyl	2	20	20.6	103		61	148