

Surrogate Summary

SW-846

SDG No.: Q3368

Client: Aramark Uniforms

Analytical Method: 8270E

Lab Sample ID	Client ID	Parameter	Spike (PPM)	Result (PPM)	Recovery (%)	Qual	Limits (%)	
							Low	High
PB170172TB	PB170172TB	2-Fluorophenol	150	146	98		23	138
		Phenol-d6	150	136	91		10	134
		Nitrobenzene-d5	100	80.4	80		67	132
		2-Fluorobiphenyl	100	77.1	77		52	132
		2,4,6-Tribromophenol	150	129	86		44	137
		Terphenyl-d14	100	89.0	89		42	152
PB170216BL	PB170216BL	2-Fluorophenol	150	151	101		23	138
		Phenol-d6	150	142	95		10	134
		Nitrobenzene-d5	100	78.7	79		67	132
		2-Fluorobiphenyl	100	76.6	77		52	132
		2,4,6-Tribromophenol	150	131	87		44	137
		Terphenyl-d14	100	80.4	80		42	152
PB170216BS	PB170216BS	2-Fluorophenol	150	140	93		23	138
		Phenol-d6	150	138	92		10	134
		Nitrobenzene-d5	100	72.9	73		67	132
		2-Fluorobiphenyl	100	69.3	69		52	132
		2,4,6-Tribromophenol	150	119	79		44	137
		Terphenyl-d14	100	83.4	83		42	152
Q3368-06	OUTSIDE-DUMPSTER	2-Fluorophenol	150	130	87		23	138
		Phenol-d6	150	122	81		10	134
		Nitrobenzene-d5	100	74.6	75		67	132
		2-Fluorobiphenyl	100	68.9	69		52	132
		2,4,6-Tribromophenol	150	129	86		44	137
		Terphenyl-d14	100	55.3	55		42	152
Q3368-06MS	OUTSIDE-DUMPSTERMS	2-Fluorophenol	150	148	99		23	138
		Phenol-d6	150	141	94		10	134
		Nitrobenzene-d5	100	85.5	85		67	132
		2-Fluorobiphenyl	100	74.5	75		52	132
		2,4,6-Tribromophenol	150	140	94		44	137
		Terphenyl-d14	100	72.8	73		42	152
Q3368-06MSD	OUTSIDE-DUMPSTERMSD	2-Fluorophenol	150	145	97		23	138
		Phenol-d6	150	137	91		10	134
		Nitrobenzene-d5	100	83.9	84		67	132
		2-Fluorobiphenyl	100	76.9	77		52	132
		2,4,6-Tribromophenol	150	125	83		44	137
		Terphenyl-d14	100	66.0	66		42	152
Q3368-07	PRESS-SLUDGE	2-Fluorophenol	150	129	86		23	138
		Phenol-d6	150	119	79		10	134
		Nitrobenzene-d5	100	75.7	76		67	132
		2-Fluorobiphenyl	100	73.5	73		52	132
		2,4,6-Tribromophenol	150	94.4	63		44	137
		Terphenyl-d14	100	58.0	58		42	152