

Surrogate Summary

SDG No.: Q3368

Client: Aramark Uniforms

Analytical Method: SW8260D

Lab Sample ID	Client ID	Parameter	Spike	Result	Recovery (%)	Qual	Limits (%)	
							Low	High
Q3368-06	OUTSIDE-DUMPSTER	1,2-Dichloroethane-d4	50	51.5	103		74	125
		Dibromofluoromethane	50	45.5	91		75	124
		Toluene-d8	50	44.8	90		86	113
		4-Bromofluorobenzene	50	50.6	101		77	121
Q3368-07	PRESS-SLUDGE	1,2-Dichloroethane-d4	50	53.8	108		74	125
		Dibromofluoromethane	50	45.1	90		75	124
		Toluene-d8	50	47.8	96		86	113
		4-Bromofluorobenzene	50	59.3	119		77	121

Surrogate Summary

SDG No.: Q3368

Client: Aramark Uniforms

Analytical Method: SW8260-Low

Lab Sample ID	Client ID	Parameter	Spike	Result	Recovery (%)	Qual	Limits (%)	
							Low	High
VX1022WBL01	VX1022WBL01	1,2-Dichloroethane-d4	50	50.7	101		74	125
		Dibromofluoromethane	50	49.1	98		75	124
		Toluene-d8	50	47.7	95		86	113
		4-Bromofluorobenzene	50	54.5	109		77	121
VX1022WBS01	VX1022WBS01	1,2-Dichloroethane-d4	50	50.2	100		74	125
		Dibromofluoromethane	50	50.9	102		75	124
		Toluene-d8	50	50.7	101		86	113
		4-Bromofluorobenzene	50	52.4	105		77	121
VX1022WBSD0	VX1022WBSD01	1,2-Dichloroethane-d4	50	49.3	99		74	125
		Dibromofluoromethane	50	49.3	99		75	124
		Toluene-d8	50	49.2	98		86	113
		4-Bromofluorobenzene	50	49.4	99		77	121