

SOP ID : M1312-SPLP-10
SDG No : N/A
Weigh By : JP
Balance ID : WC SC-7
pH Meter ID : WC PH METER-1
Extraction By : N/A
Filter By : JP
Pipette ID : WC
Tumbler ID : T-1
TCLP Filter ID : 40718710

Start Prep Date : 10/28/2025 **Time :** 15:00
End Prep Date : 10/29/2025 **Time :** 08:35
Combination Ratio : 20
ZHE Cleaning Batch : N/A 10
Initial Room Temperature: 24 °C
Final Room Temperature: 22 °C
TCLP Technician Signature : 10
Supervisor By : 12

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
SPLP FLUID	WP114530	N/A
N/A	N/A	N/A
N/A	N/A	N/A
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1940,W1941,W1942	W3166,W1938,W1939,
N/A	N/A	N/A
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP86978

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. q3391-01 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
10/29/25 10:00	10 1000m	5125 1000m
	Preparation Group	Analysis Group

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB170309TB	LEB309	02	N/A	2000	N/A	N/A	N/A	4.24	1.5	T-1
Q3391-01	RCA	01	100.03	2000	N/A	N/A	N/A	12.0	1.0	T-1

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB170309TB	LEB309	N/A	N/A	N/A	N/A	N/A	N/A
Q3391-01	RCA	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 / WC S-2

Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB170309TB	LEB309	N/A	N/A	N/A	N/A	#1	4.24
Q3391-01	RCA	N/A	N/A	N/A	N/A	#1	4.24