

### Matrix Spike/Matrix Spike Duplicate Summary

SW-846

|                 |                                |                           |                   |
|-----------------|--------------------------------|---------------------------|-------------------|
| <b>SDG No.:</b> | <u>Q3399</u>                   | <b>Analytical Method:</b> | <u>8151A</u>      |
| <b>Client:</b>  | <u>ICE Service Group, Inc.</u> | <b>DataFile :</b>         | <u>PS032145.D</u> |

|                       |                   | Sample                   |        |                       |       | Rec |      | RPD |      | Limits |      |
|-----------------------|-------------------|--------------------------|--------|-----------------------|-------|-----|------|-----|------|--------|------|
|                       | Parameter         | Spike                    | Result | Result                | Units | Rec | Qual | RPD | Qual | Low    | High |
| <b>Lab Sample ID:</b> | <b>Q3399-02MS</b> | <b>Client Sample ID:</b> |        | <b>COMP-ROLLOFFMS</b> |       |     |      |     |      |        |      |
|                       | (Column 1)        |                          |        |                       |       |     |      |     |      |        |      |
|                       | 2,4-D             | 50                       | 0      | 25.1                  | ug/L  | 50  |      |     |      | 45     | 152  |
|                       | 2,4,5-TP(Silvex)  | 50                       | 0      | 24.5                  | ug/L  | 49  | *    |     |      | 51     | 134  |
| <b>Lab Sample ID:</b> | <b>Q3399-02MS</b> | <b>Client Sample ID:</b> |        | <b>COMP-ROLLOFFMS</b> |       |     |      |     |      |        |      |
|                       | (Column 2)        |                          |        |                       |       |     |      |     |      |        |      |
|                       | 2,4-D             | 50                       | 0      | 29.1                  | ug/L  | 58  |      |     |      | 45     | 152  |
|                       | 2,4,5-TP(Silvex)  | 50                       | 0      | 30.6                  | ug/L  | 61  |      |     |      | 51     | 134  |