

Instrument ID: BNA_F

Daily Analysis Runlog For Sequence/QC Batch ID # BF102425

Review By	rahul	Review On	10/27/2025 10:16:31 AM
Supervise By	Jagrut	Supervise On	10/27/2025 4:26:19 PM
SubDirectory	BF102425	HP Acquire Method	BNA_F
		HP Processing Method	BF102425

STD. NAME	STD REF.#
Tune/Reschk	SP6875
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864
CCC	SP6860
Internal Standard/PEM	S13179,10ul/1000ul sample
ICV/I.BLK	SP6872
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BF144055.D	24 Oct 2025 09:45		RC/JU	Ok
2	SSTDICC2.5	SSTDICC2.5	BF144056.D	24 Oct 2025 10:13		RC/JU	Ok
3	SSTDICC005	SSTDICC005	BF144057.D	24 Oct 2025 10:42		RC/JU	Ok
4	SSTDICC010	SSTDICC010	BF144058.D	24 Oct 2025 11:39		RC/JU	Ok
5	SSTDICC020	SSTDICC020	BF144059.D	24 Oct 2025 12:36		RC/JU	Ok
6	SSTDICCC040	SSTDICCC040	BF144060.D	24 Oct 2025 13:05		RC/JU	Ok
7	SSTDICC050	SSTDICC050	BF144061.D	24 Oct 2025 13:34		RC/JU	Ok
8	SSTDICC060	SSTDICC060	BF144062.D	24 Oct 2025 14:03	Not Used	RC/JU	Not Ok
9	SSTDICC080	SSTDICC080	BF144063.D	24 Oct 2025 14:32		RC/JU	Ok
10	SSTDICC060	SSTDICC060	BF144064.D	24 Oct 2025 16:24		RC/JU	Ok
11	SSTDICV040	ICVBF102425	BF144065.D	24 Oct 2025 17:22		RC/JU	Ok
12	PB170222BL	PB170222BL	BF144066.D	24 Oct 2025 17:51		RC/JU	Ok
13	Q3431-01	EO-02-10222025	BF144067.D	24 Oct 2025 18:20		RC/JU	Ok,M
14	Q3431-01DL	EO-02-10222025DL	BF144068.D	24 Oct 2025 18:48	Not Required	RC/JU	Not Ok
15	Q3405-01	PAD-10-20-25-A	BF144069.D	24 Oct 2025 19:17		RC/JU	Ok
16	Q3406-01	ARS20-0037	BF144070.D	24 Oct 2025 19:46		RC/JU	Ok
17	Q3412-01	TP-6	BF144071.D	24 Oct 2025 20:15	Need 5X Dilution	RC/JU	Dilution
18	Q3445-01	EO-01-102325	BF144072.D	24 Oct 2025 20:44		RC/JU	Ok,M

Instrument ID: BNA_F

Daily Analysis Runlog For Sequence/QC Batch ID # BF102425

Review By	rahul	Review On	10/27/2025 10:16:31 AM		
Supervise By	Jagrut	Supervise On	10/27/2025 4:26:19 PM		
SubDirectory	BF102425	HP Acquire Method	BNA_F	HP Processing Method	BF102425
STD. NAME	STD REF.#				
Tune/Reschk	SP6875				
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864				
CCC	SP6860				
Internal Standard/PEM	S13179,10ul/1000ul sample				
ICV/I.BLK	SP6872				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration