

Instrument ID: BNA_F

Daily Analysis Runlog For Sequence/QC Batch ID # BF102825

Review By	rahul	Review On	10/28/2025 3:29:21 PM
Supervise By	Jagrut	Supervise On	10/29/2025 3:24:40 PM
SubDirectory	BF102825	HP Acquire Method	BNA_F
		HP Processing Method	BF102425

STD. NAME	STD REF.#
Tune/Reschk	SP6875
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864
CCC	SP6860
Internal Standard/PEM	S13180,10ul/1000ul sample
ICV/I.BLK	SP6872
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BF144096.D	28 Oct 2025 12:48		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BF144097.D	28 Oct 2025 13:17		RC/JU	Ok
3	PB170267BL	PB170267BL	BF144098.D	28 Oct 2025 13:47		RC/JU	Ok
4	PB170267BS	PB170267BS	BF144099.D	28 Oct 2025 14:16		RC/JU	Ok,M
5	PB170202TB	PB170202TB	BF144100.D	28 Oct 2025 14:46	Internal standard fail	RC/JU	Not Ok
6	Q3424-01	PR132-DP3-20251021	BF144101.D	28 Oct 2025 15:16		RC/JU	Ok
7	Q3424-02	PR132-S11-010080-20	BF144102.D	28 Oct 2025 15:45		RC/JU	Ok
8	Q3424-03	PR132-S14-000094-20	BF144103.D	28 Oct 2025 16:15		RC/JU	Ok
9	Q3395-01	PR132-S09-000051-20	BF144104.D	28 Oct 2025 16:44		RC/JU	Ok
10	Q3419-04	PR132-S16-119133-20	BF144105.D	28 Oct 2025 17:14	Hit for pyrene & it is biased high in the CCC.	RC/JU	ReRun
11	Q3434-01	PR132-WC1-20251022	BF144106.D	28 Oct 2025 17:44	Hit for pyrene below CRQL	RC/JU	Ok
12	Q3460-09	WC-3	BF144107.D	28 Oct 2025 18:14		RC/JU	Ok
13	Q3439-08	TP-10	BF144108.D	28 Oct 2025 18:44		RC/JU	Ok
14	Q3440-04	JB-2	BF144109.D	28 Oct 2025 19:13		RC/JU	Ok
15	Q3460-04	WC-1	BF144110.D	28 Oct 2025 19:43		RC/JU	Ok
16	Q3460-08	WC-2	BF144111.D	28 Oct 2025 20:13		RC/JU	Ok
17	Q3460-12	WC-3	BF144112.D	28 Oct 2025 20:43		RC/JU	Ok

Instrument ID: BNA_F

Daily Analysis Runlog For Sequence/QC Batch ID # BF102825

Review By	rahul	Review On	10/28/2025 3:29:21 PM		
Supervise By	Jagrut	Supervise On	10/29/2025 3:24:40 PM		
SubDirectory	BF102825	HP Acquire Method	BNA_F	HP Processing Method	BF102425
STD. NAME	STD REF.#				
Tune/Reschk	SP6875				
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864				
CCC	SP6860				
Internal Standard/PEM	S13180,10ul/1000ul sample				
ICV/I.BLK	SP6872				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

18	Q3427-04	TP-8	BF144113.D	28 Oct 2025 21:12		RC/JU	Ok
19	Q3439-04	TP-9	BF144114.D	28 Oct 2025 21:42		RC/JU	Ok
20	Q3395-06	PR132-COMP03-2025	BF144115.D	28 Oct 2025 22:11		RC/JU	Ok
21	Q3447-01DL	P001-TW1-251023-01	BF144116.D	28 Oct 2025 22:40		RC/JU	Ok
22	Q3395-05	PR132-S07-102116-20	BF144117.D	28 Oct 2025 23:10		RC/JU	Ok
23	Q3451-08	TP-12	BF144118.D	28 Oct 2025 23:39		RC/JU	Ok
24	Q3447-06	P001-TW3-251023-01	BF144119.D	29 Oct 2025 00:09		RC/JU	Ok

M : Manual Integration