

Instrument ID: BNA_G

Daily Analysis Runlog For Sequence/QC Batch ID # BG102825

Review By	rahul	Review On	10/28/2025 3:29:34 PM
Supervise By	mohammad	Supervise On	11/5/2025 12:11:31 AM
SubDirectory	BG102825	HP Acquire Method	BNA_G
		HP Processing Method	BG102425

STD. NAME	STD REF.#
Tune/Reschk	SP6875
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864
CCC	SP6860
Internal Standard/PEM	S13180,10ul/1000ul sample
ICV/I.BLK	SP6872
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BG064599.D	28 Oct 2025 13:05		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BG064600.D	28 Oct 2025 13:45		RC/JU	Ok
3	PB170239TB	PB170239TB	BG064601.D	28 Oct 2025 14:26		RC/JU	Ok
4	Q3460-16	WC-4	BG064602.D	28 Oct 2025 15:10		RC/JU	Ok
5	Q3460-20	WC-5	BG064603.D	28 Oct 2025 15:51		RC/JU	Ok
6	Q3452-04	JB-1	BG064604.D	28 Oct 2025 16:32		RC/JU	Ok
7	Q3451-04	TP-11	BG064605.D	28 Oct 2025 17:13		RC/JU	Ok
8	Q3419-05	PR132-S16-094119-20	BG064606.D	28 Oct 2025 17:54		RC/JU	Ok,M
9	Q3427-01	TP-8	BG064607.D	28 Oct 2025 18:35		RC/JU	Ok,M
10	Q3452-01	JB-1	BG064608.D	28 Oct 2025 19:16		RC/JU	Ok,M
11	Q3376-05	PR132-S12-020047-20	BG064609.D	28 Oct 2025 19:57		RC/JU	Ok
12	Q3451-05	TP-12	BG064610.D	28 Oct 2025 20:37		RC/JU	Ok,M
13	Q3460-05	WC-2	BG064611.D	28 Oct 2025 21:18		RC/JU	Ok
14	Q3460-01	WC-1	BG064612.D	28 Oct 2025 21:59		RC/JU	Ok,M
15	Q3460-13	WC-4	BG064613.D	28 Oct 2025 22:39		RC/JU	Ok,M
16	Q3460-17	WC-5	BG064614.D	28 Oct 2025 23:20	Need Straight Run	RC/JU	Not Ok
17	Q3457-01	TR-05-102425	BG064615.D	29 Oct 2025 00:01		RC/JU	Ok,M
18	Q3465-01	SU-03-10242025	BG064616.D	29 Oct 2025 00:41		RC/JU	Ok

Instrument ID: BNA_G

Daily Analysis Runlog For Sequence/QC Batch ID # BG102825

Review By	rahul	Review On	10/28/2025 3:29:34 PM		
Supervise By	mohammad	Supervise On	11/5/2025 12:11:31 AM		
SubDirectory	BG102825	HP Acquire Method	BNA_G	HP Processing Method	BG102425
STD. NAME	STD REF.#				
Tune/Reschk	SP6875				
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864				
CCC	SP6860				
Internal Standard/PEM	S13180,10ul/1000ul sample				
ICV/I.BLK	SP6872				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration