

### Matrix Spike/Matrix Spike Duplicate Summary

SW-846

|                 |                                     |                           |                   |
|-----------------|-------------------------------------|---------------------------|-------------------|
| <b>SDG No.:</b> | <u>Q3437</u>                        | <b>Analytical Method:</b> | <u>8082A</u>      |
| <b>Client:</b>  | <u>Matrix New World Engineering</u> | <b>DataFile :</b>         | <u>PP076281.D</u> |

|                | Parameter                 | Spike             | Sample |                     | Units | Rec | Rec  |     | RPD  |         | Limits    |         |
|----------------|---------------------------|-------------------|--------|---------------------|-------|-----|------|-----|------|---------|-----------|---------|
|                |                           |                   | Result | Result              |       |     | Qual | RPD | Qual | Low     | High      | RPD     |
| Lab Sample ID: | Q3450-04MSD<br>(Column 1) | Client Sample ID: |        | SS-10R-(3.5-4.0)MSD |       |     |      |     |      |         |           |         |
|                | AR1016                    | 192.2             | 0      | 189                 | ug/kg | 98  |      | 2   |      | 40 (55) | 140 (146) | 30 (15) |
|                | AR1260                    | 192.2             | 0      | 215                 | ug/kg | 112 |      | 3   |      | 40 (54) | 140 (119) | 30 (15) |
| Lab Sample ID: | Q3450-04MSD<br>(Column 2) | Client Sample ID: |        | SS-10R-(3.5-4.0)MSD |       |     |      |     |      |         |           |         |
|                | AR1016                    | 192.2             | 0      | 173                 | ug/kg | 90  |      | 1   |      | 40 (55) | 140 (146) | 30 (15) |
|                | AR1260                    | 192.2             | 0      | 206                 | ug/kg | 107 |      | 2   |      | 40 (54) | 140 (119) | 30 (15) |