

Surrogate Summary

SDG No.: Q3447

Client: Weston Solutions, Inc.

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PP075887.D	PIBLK-PP075887.D	Tetrachloro-m-xyl	1	20	19.9	100		60	140
		Decachlorobiphen	1	20	20.4	102		60	140
		Tetrachloro-m-xyl	2	20	18.8	94		60	140
		Decachlorobiphen	2	20	19.3	96		60	140
I.BLK-PP076067.D	PIBLK-PP076067.D	Tetrachloro-m-xyl	1	20	22.1	110		60	140
		Decachlorobiphen	1	20	22.4	112		60	140
		Tetrachloro-m-xyl	2	20	18.9	95		60	140
		Decachlorobiphen	2	20	18.2	91		60	140
PB170236BL	PB170236BL	Tetrachloro-m-xyl	1	20	27.3	137		32	144
		Decachlorobiphen	1	20	29.0	145		32	175
		Tetrachloro-m-xyl	2	20	26.1	130		32	144
		Decachlorobiphen	2	20	25.2	126		32	175
PB170236BS	PB170236BS	Tetrachloro-m-xyl	1	20	26.7	134		32	144
		Decachlorobiphen	1	20	28.6	143		32	175
		Tetrachloro-m-xyl	2	20	21.3	107		32	144
		Decachlorobiphen	2	20	24.3	121		32	175
Q3447-03	P001-TW2-251023-01	Tetrachloro-m-xyl	1	20	40.9	205	*	32	144
		Decachlorobiphen	1	20	32.9	165		32	175
		Tetrachloro-m-xyl	2	20	265	1324	*	32	144
		Decachlorobiphen	2	20	26.1	131		32	175
Q3447-07	P001-TW4-251023-01	Tetrachloro-m-xyl	1	20	11.0	55		32	144
		Decachlorobiphen	1	20	34.2	171		32	175
		Tetrachloro-m-xyl	2	20	92.7	463	*	32	144
		Decachlorobiphen	2	20	7.20	36		32	175
I.BLK-PP076082.D	PIBLK-PP076082.D	Tetrachloro-m-xyl	1	20	18.9	95		60	140
		Decachlorobiphen	1	20	20.4	102		60	140
		Tetrachloro-m-xyl	2	20	17.1	86		60	140
		Decachlorobiphen	2	20	16.7	83		60	140
Q3447-05	P001-TW3-251023-01	Tetrachloro-m-xyl	1	20	27.3	137		32	144
		Decachlorobiphen	1	20	30.0	150		32	175
		Tetrachloro-m-xyl	2	20	33.3	167	*	32	144
		Decachlorobiphen	2	20	26.6	133		32	175
Q3447-01	P001-TW1-251023-01	Tetrachloro-m-xyl	1	20	115	575	*	32	144
		Decachlorobiphen	1	20	26.5	133		32	175
		Tetrachloro-m-xyl	2	20	934	4668	*	32	144
		Decachlorobiphen	2	20	12.0	60		32	175
I.BLK-PP076097.D	PIBLK-PP076097.D	Tetrachloro-m-xyl	1	20	18.7	93		60	140
		Decachlorobiphen	1	20	20.4	102		60	140
		Tetrachloro-m-xyl	2	20	16.8	84		60	140
		Decachlorobiphen	2	20	17.5	87		60	140
I.BLK-PP076115.D	PIBLK-PP076115.D	Tetrachloro-m-xyl	1	20	18.3	92		60	140

Surrogate Summary

SDG No.: Q3447

Client: Weston Solutions, Inc.

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PP076115.D	PIBLK-PP076115.D	Decachlorobiphen	1	20	19.4	97		60	140
		Tetrachloro-m-xyl	2	20	16.2	81		60	140
Q3447-03RE	P001-TW2-251023-01RE	Decachlorobiphen	2	20	15.9	79		60	140
		Tetrachloro-m-xyl	1	20	44.3	221	*	32	144
		Decachlorobiphen	1	20	33.8	169		32	175
		Tetrachloro-m-xyl	2	20	264	1318	*	32	144
		Decachlorobiphen	2	20	23.7	119		32	175
Q3447-01RE	P001-TW1-251023-01RE	Tetrachloro-m-xyl	1	20	134	670	*	32	144
		Decachlorobiphen	1	20	141	705	*	32	175
		Tetrachloro-m-xyl	2	20	931	4655	*	32	144
		Decachlorobiphen	2	20	11.0	55		32	175
I.BLK-PP076131.D	PIBLK-PP076131.D	Tetrachloro-m-xyl	1	20	17.9	90		60	140
		Decachlorobiphen	1	20	18.8	94		60	140
		Tetrachloro-m-xyl	2	20	16.3	81		60	140
		Decachlorobiphen	2	20	16.2	81		60	140
I.BLK-PQ070972.D	PIBLK-PQ070972.D	Tetrachloro-m-xyl	1	20	18.3	91		60	140
		Decachlorobiphen	1	20	19.0	95		60	140
		Tetrachloro-m-xyl	2	20	19.1	96		60	140
		Decachlorobiphen	2	20	19.4	97		60	140
I.BLK-PQ071150.D	PIBLK-PQ071150.D	Tetrachloro-m-xyl	1	20	22.5	113		60	140
		Decachlorobiphen	1	20	25.7	129		60	140
		Tetrachloro-m-xyl	2	20	19.9	100		60	140
		Decachlorobiphen	2	20	24.5	122		60	140
Q3450-02MS	SS-5R-(3.5-4.0)MS	Tetrachloro-m-xyl	1	20	23.4	117		32	144
		Decachlorobiphen	1	20	22.4	112		32	175
		Tetrachloro-m-xyl	2	20	18.7	94		32	144
		Decachlorobiphen	2	20	20.7	104		32	175
Q3450-02MSD	SS-5R-(3.5-4.0)MSD	Tetrachloro-m-xyl	1	20	25.7	128		32	144
		Decachlorobiphen	1	20	24.0	120		32	175
		Tetrachloro-m-xyl	2	20	22.2	111		32	144
		Decachlorobiphen	2	20	22.9	114		32	175
I.BLK-PQ071165.D	PIBLK-PQ071165.D	Tetrachloro-m-xyl	1	20	20.3	102		60	140
		Decachlorobiphen	1	20	23.5	117		60	140
		Tetrachloro-m-xyl	2	20	18.4	92		60	140
		Decachlorobiphen	2	20	22.2	111		60	140