

Instrument ID: BNA_G

Daily Analysis Runlog For Sequence/QC Batch ID # BG102925

Review By	rahul	Review On	10/30/2025 8:51:46 AM
Supervise By	Jagrut	Supervise On	10/30/2025 4:39:08 PM
SubDirectory	BG102925	HP Acquire Method	BNA_G
		HP Processing Method	BG102425

STD. NAME	STD REF.#
Tune/Reschk	SP6875
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864
CCC	SP6860
Internal Standard/PEM	S13180,10ul/1000ul sample
ICV/I.BLK	SP6872
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BG064617.D	29 Oct 2025 10:41		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BG064618.D	29 Oct 2025 11:22		RC/JU	Ok
3	PB170285BL	PB170285BL	BG064619.D	29 Oct 2025 12:03		RC/JU	Ok
4	PB170285BS	PB170285BS	BG064620.D	29 Oct 2025 12:44		RC/JU	Ok,M
5	PB170215BL	PB170215BL	BG064621.D	29 Oct 2025 13:25		RC/JU	Ok
6	PB170215BS	PB170215BS	BG064622.D	29 Oct 2025 14:06		RC/JU	Ok,M
7	Q3477-01	SOUTH-MAHWAH-SOI	BG064623.D	29 Oct 2025 14:50		RC/JU	Ok
8	Q3477-01MS	SOUTH-MAHWAH-SOI	BG064624.D	29 Oct 2025 15:31		RC/JU	Ok,M
9	Q3477-01MSD	SOUTH-MAHWAH-SOI	BG064625.D	29 Oct 2025 16:12		RC/JU	Ok,M
10	Q3419-01	PR-132-S16-015094-20	BG064626.D	29 Oct 2025 16:53		RC/JU	Ok
11	Q3419-02MS	PR-132-S16-015094-20	BG064627.D	29 Oct 2025 17:34		RC/JU	Ok,M
12	Q3419-03MSD	PR-132-S16-015094-20	BG064628.D	29 Oct 2025 18:15		RC/JU	Ok,M
13	Q3468-01	MW-06	BG064629.D	29 Oct 2025 18:56		RC/JU	Ok,M
14	Q3468-02	MW-08	BG064630.D	29 Oct 2025 19:37		RC/JU	Ok,M
15	Q3468-03	MW-10	BG064631.D	29 Oct 2025 20:18		RC/JU	Ok
16	Q3468-04	MW-11	BG064632.D	29 Oct 2025 20:59		RC/JU	Ok
17	Q3468-05	MW-13	BG064633.D	29 Oct 2025 21:40		RC/JU	Ok
18	Q3468-06	MW-12	BG064634.D	29 Oct 2025 22:21		RC/JU	Ok

Instrument ID: BNA_G

Daily Analysis Runlog For Sequence/QC Batch ID # BG102925

Review By	rahul	Review On	10/30/2025 8:51:46 AM		
Supervise By	Jagrut	Supervise On	10/30/2025 4:39:08 PM		
SubDirectory	BG102925	HP Acquire Method	BNA_G	HP Processing Method	BG102425
STD. NAME	STD REF.#				
Tune/Reschk	SP6875				
Initial Calibration Stds	SP6857,SP6858,SP6859,SP6860,SP6861,SP6862,SP6863,SP6864				
CCC	SP6860				
Internal Standard/PEM	S13180,10ul/1000ul sample				
ICV/I.BLK	SP6872				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration