

Surrogate Summary

SDG No.: Q3495

Client: Remington & Vernick Engineers

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO114312.D	PIBLK-PO114312.D	Tetrachloro-m-xyl	1	20	18.6	93		60	140
		Decachlorobiphen	1	20	19.2	96		60	140
		Tetrachloro-m-xyl	2	20	18.9	94		60	140
		Decachlorobiphen	2	20	19.0	95		60	140
I.BLK-PO114802.D	PIBLK-PO114802.D	Tetrachloro-m-xyl	1	20	18.5	92		60	140
		Decachlorobiphen	1	20	17.9	90		60	140
		Tetrachloro-m-xyl	2	20	16.4	82		60	140
		Decachlorobiphen	2	20	19.5	98		60	140
PB170339BL	PB170339BL	Tetrachloro-m-xyl	1	20	27.3	136		32	144
		Decachlorobiphen	1	20	26.8	134		32	175
		Tetrachloro-m-xyl	2	20	24.8	124		32	144
		Decachlorobiphen	2	20	26.2	131		32	175
PB170339BS	PB170339BS	Tetrachloro-m-xyl	1	20	27.4	137		32	144
		Decachlorobiphen	1	20	26.7	134		32	175
		Tetrachloro-m-xyl	2	20	23.1	115		32	144
		Decachlorobiphen	2	20	26.3	132		32	175
Q3495-02	SB-1025-2	Tetrachloro-m-xyl	1	20	28.2	141		32	144
		Decachlorobiphen	1	20	21.7	109		32	175
		Tetrachloro-m-xyl	2	20	25.5	128		32	144
		Decachlorobiphen	2	20	22.3	112		32	175
Q3495-02MS	SB-1025-2MS	Tetrachloro-m-xyl	1	20	21.6	108		32	144
		Decachlorobiphen	1	20	15.3	76		32	175
		Tetrachloro-m-xyl	2	20	18.4	92		32	144
		Decachlorobiphen	2	20	16.4	82		32	175
Q3495-02MSD	SB-1025-2MSD	Tetrachloro-m-xyl	1	20	21.7	109		32	144
		Decachlorobiphen	1	20	15.8	79		32	175
		Tetrachloro-m-xyl	2	20	18.6	93		32	144
		Decachlorobiphen	2	20	16.9	84		32	175
I.BLK-PO114817.D	PIBLK-PO114817.D	Tetrachloro-m-xyl	1	20	17.5	87		60	140
		Decachlorobiphen	1	20	17.6	88		60	140
		Tetrachloro-m-xyl	2	20	16.0	80		60	140
		Decachlorobiphen	2	20	18.8	94		60	140
I.BLK-PP075887.D	PIBLK-PP075887.D	Tetrachloro-m-xyl	1	20	19.9	100		60	140
		Decachlorobiphen	1	20	20.4	102		60	140
		Tetrachloro-m-xyl	2	20	18.8	94		60	140
		Decachlorobiphen	2	20	19.3	96		60	140
I.BLK-PP076260.D	PIBLK-PP076260.D	Tetrachloro-m-xyl	1	20	22.0	110		60	140
		Decachlorobiphen	1	20	21.4	107		60	140
		Tetrachloro-m-xyl	2	20	19.2	96		60	140
		Decachlorobiphen	2	20	17.9	90		60	140
PB170347BL	PB170347BL	Tetrachloro-m-xyl	1	20	32.5	162		30	173

Surrogate Summary

SDG No.: Q3495

Client: Remington & Vernick Engineers

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
PB170347BL	PB170347BL	Decachlorobiphen	1	20	31.3	157		10	173
		Tetrachloro-m-xyl	2	20	29.0	145		30	173
PB170347BS	PB170347BS	Decachlorobiphen	2	20	26.5	133		10	173
		Tetrachloro-m-xyl	1	20	32.9	165		30	173
		Decachlorobiphen	1	20	30.7	153		10	173
		Tetrachloro-m-xyl	2	20	29.5	147		30	173
		Decachlorobiphen	2	20	26.0	130		10	173
PB170347BSD	PB170347BSD	Tetrachloro-m-xyl	1	20	31.0	155		30	173
		Decachlorobiphen	1	20	29.8	149		10	173
		Tetrachloro-m-xyl	2	20	29.4	147		30	173
		Decachlorobiphen	2	20	25.0	125		10	173
Q3495-03	TW-1025-1	Tetrachloro-m-xyl	1	20	27.2	136		30	173
		Decachlorobiphen	1	20	21.3	106		10	173
		Tetrachloro-m-xyl	2	20	21.0	105		30	173
		Decachlorobiphen	2	20	18.1	91		10	173
I.BLK-PP076275.D	PIBLK-PP076275.D	Tetrachloro-m-xyl	1	20	21.6	108		60	140
		Decachlorobiphen	1	20	21.0	105		60	140
		Tetrachloro-m-xyl	2	20	19.3	97		60	140
		Decachlorobiphen	2	20	17.5	87		60	140