

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE112725AL

Review By	rahul	Review On	11/26/2025 3:46:33 PM
Supervise By	mohammad	Supervise On	12/2/2025 12:49:44 AM
SubDirectory	FE112725AL	HP Acquire Method	HP Processing Method FE112725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24768,PP24770,PP24771,PP24772,PP24773		
CCC Internal Standard/PEM	PP24771		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24769,PP24774		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FE057077.D	27 Nov 2025 11:16	YP\AJ	Ok
2	I.BLK	FE057078.D	27 Nov 2025 11:47	YP\AJ	Ok
3	100 PPM ALIPHATIC HC STD1	FE057079.D	27 Nov 2025 12:17	YP\AJ	Ok
4	50 PPM ALIPHATIC HC STD2	FE057080.D	27 Nov 2025 12:47	YP\AJ	Ok
5	20 PPM ALIPHATIC HC STD3	FE057081.D	27 Nov 2025 13:17	YP\AJ	Ok
6	10 PPM ALIPHATIC HC STD4	FE057082.D	27 Nov 2025 13:47	YP\AJ	Ok,M
7	5 PPM ALIPHATIC HC STD5	FE057083.D	27 Nov 2025 14:17	YP\AJ	Ok,M
8	20 PPM ALIPHATIC HC STD ICV	FE057084.D	27 Nov 2025 14:47	YP\AJ	Ok
9	I.BLK	FE057085.D	27 Nov 2025 15:48	YP\AJ	Ok
10	20 PPM ALIPHATIC HC STD	FE057086.D	27 Nov 2025 16:18	YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE121125AL

Review By	yogesh	Review On	12/11/2025 11:55:13 AM
Supervise By	rahul	Supervise On	12/12/2025 10:05:18 AM
SubDirectory	FE121125AL	HP Acquire Method	HP Processing Method FE112725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24768,PP24770,PP24771,PP24772,PP24773		
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24771 PP24769,PP24774		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FE057306.D	11 Dec 2025 08:24	YP\AJ	Ok
2	I.BLK	FE057307.D	11 Dec 2025 08:55	YP\AJ	Ok
3	20 PPM ALIPHATIC HC STD	FE057308.D	11 Dec 2025 09:25	YP\AJ	Ok
4	PB170903BL	FE057309.D	11 Dec 2025 13:32	YP\AJ	Ok
5	PB170903BS	FE057310.D	11 Dec 2025 14:02	YP\AJ	Ok,M
6	PB170903BSD	FE057311.D	11 Dec 2025 14:33	YP\AJ	Ok,M
7	Q3831-01	FE057312.D	11 Dec 2025 15:03	YP\AJ	Ok
8	Q3831-01D	FE057313.D	11 Dec 2025 15:33	YP\AJ	Ok
9	Q3831-01MS	FE057314.D	11 Dec 2025 16:03	YP\AJ	Ok,M
10	Q3831-01MSD	FE057315.D	11 Dec 2025 16:34	YP\AJ	Ok,M
11	Q3831-02	FE057316.D	11 Dec 2025 17:04	YP\AJ	Ok
12	Q3831-03	FE057317.D	11 Dec 2025 17:34	YP\AJ	Ok
13	I.BLK	FE057318.D	11 Dec 2025 18:35	YP\AJ	Ok
14	20 PPM ALIPHATIC HC STD	FE057319.D	11 Dec 2025 19:05	YP\AJ	Ok
15	PB170799BL	FE057320.D	11 Dec 2025 20:06	YP\AJ	Ok
16	PB170800BL	FE057321.D	11 Dec 2025 20:36	YP\AJ	Ok
17	Q3530-01	FE057322.D	11 Dec 2025 21:06	YP\AJ	Ok
18	Q3530-01	FE057323.D	11 Dec 2025 21:37	YP\AJ	Ok
19	Q3530-02	FE057324.D	11 Dec 2025 22:07	YP\AJ	Ok
20	Q3530-07	FE057325.D	11 Dec 2025 22:37	YP\AJ	Ok
21	Q3530-07	FE057326.D	11 Dec 2025 23:07	YP\AJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE121125AL

Review By	yogesh	Review On	12/11/2025 11:55:13 AM
Supervise By	rahul	Supervise On	12/12/2025 10:05:18 AM
SubDirectory	FE121125AL	HP Acquire Method	HP Processing Method FE112725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24768,PP24770,PP24771,PP24772,PP24773		
CCC Internal Standard/PEM	PP24771		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24769,PP24774		

22	Q3530-08	FE057327.D	11 Dec 2025 23:37	YP\AJ	Ok
23	I.BLK	FE057328.D	12 Dec 2025 00:38	YP\AJ	Ok
24	20 PPM ALIPHATIC HC STD	FE057329.D	12 Dec 2025 01:08	YP\AJ	Ok

M : Manual Integration