

**metals**  
**- 14 -**  
**ANALYSIS RUN LOG**

**Client:** Spectra East Inc.

**Contract:** SPEC01

**Lab code:** ACE

**Sdg no.:** Q3551

**Instrument id number:** \_\_\_\_\_

**Method:** \_\_\_\_\_

**Run number:** LB137819

**Start date:** 11/07/2025

**End date:** 11/07/2025

| Lab sample id. | Client Sample Id | d/f | Time | Parameter list |
|----------------|------------------|-----|------|----------------|
| S0             | S0               | 1   | 1351 | HG             |
| S0.05          | S0.05            | 1   | 1353 | HG             |
| S0.2           | S0.2             | 1   | 1355 | HG             |
| S2.5           | S2.5             | 1   | 1358 | HG             |
| S5             | S5               | 1   | 1400 | HG             |
| S7.5           | S7.5             | 1   | 1417 | HG             |
| S10            | S10              | 1   | 1423 | HG             |
| ICV103         | ICV103           | 1   | 1428 | HG             |
| ICB103         | ICB103           | 1   | 1431 | HG             |
| CCV40          | CCV40            | 1   | 1433 | HG             |
| CCB40          | CCB40            | 1   | 1435 | HG             |
| CRA            | CRA              | 1   | 1443 | HG             |
| PB170437BL     | PB170437BL       | 1   | 1450 | HG             |
| PB170437BS     | PB170437BS       | 1   | 1452 | HG             |
| Q3551-01       | MANHOLE          | 1   | 1457 | HG             |
| Q3551-01DUP    | MANHOLEDUP       | 1   | 1459 | HG             |
| Q3551-01MS     | MANHOLEMS        | 1   | 1501 | HG             |
| Q3551-01MSD    | MANHOLEMSD       | 1   | 1504 | HG             |
| CCV41          | CCV41            | 1   | 1506 | HG             |
| CCB41          | CCB41            | 1   | 1508 | HG             |
| Q3551-01L      | MANHOLEL         | 5   | 1511 | HG             |
| CCV42          | CCV42            | 1   | 1515 | HG             |
| CCB42          | CCB42            | 1   | 1517 | HG             |

**metals**  
**- 14 -**  
**ANALYSIS RUN LOG**

**Client:** Spectra East Inc.

**Contract:** SPEC01

**Lab code:** ACE

**Sdg no.:** Q3551

**Instrument id number:** \_\_\_\_\_

**Method:** \_\_\_\_\_

**Run number:** LB137841

**Start date:** 11/10/2025

**End date:** 11/10/2025

| Lab sample id. | Client Sample Id | d/f | Time | Parameter list                |
|----------------|------------------|-----|------|-------------------------------|
| S0             | S0               | 1   | 1120 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| S1             | S1               | 1   | 1124 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| S2             | S2               | 1   | 1128 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| S3             | S3               | 1   | 1132 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| S4             | S4               | 1   | 1137 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| S5             | S5               | 1   | 1141 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICV01          | ICV01            | 1   | 1145 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| LLICV01        | LLICV01          | 1   | 1157 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICB01          | ICB01            | 1   | 1211 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CRI01          | CRI01            | 1   | 1216 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICSA01         | ICSA01           | 1   | 1221 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICSAB01        | ICSAB01          | 1   | 1226 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICSA           | ICSA             | 20  | 1230 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| ICSAB          | ICSAB            | 20  | 1234 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV01          | CCV01            | 1   | 1238 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB01          | CCB01            | 1   | 1242 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV02          | CCV02            | 1   | 1328 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB02          | CCB02            | 1   | 1354 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| PB170412BL     | PB170412BL       | 1   | 1358 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| PB170412BS     | PB170412BS       | 1   | 1403 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| Q3547-02DUP    | CompDUP          | 1   | 1417 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| Q3547-02L      | CompL            | 5   | 1421 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| Q3547-02MS     | CompMS           | 1   | 1426 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| Q3547-02MSD    | CompMSD          | 1   | 1430 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| Q3547-02A      | CompA            | 1   | 1434 | Ag                            |
| Q3551-01       | MANHOLE          | 1   | 1438 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV03          | CCV03            | 1   | 1446 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB03          | CCB03            | 1   | 1450 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV04          | CCV04            | 1   | 1544 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB04          | CCB04            | 1   | 1547 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV05          | CCV05            | 1   | 1635 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB05          | CCB05            | 1   | 1638 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV06          | CCV06            | 1   | 1727 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB06          | CCB06            | 1   | 1731 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV07          | CCV07            | 1   | 1914 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB07          | CCB07            | 1   | 1918 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV08          | CCV08            | 1   | 2011 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB08          | CCB08            | 1   | 2026 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV09          | CCV09            | 1   | 2131 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB09          | CCB09            | 1   | 2141 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV10          | CCV10            | 1   | 2229 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |

**metals**  
**- 14 -**  
**ANALYSIS RUN LOG**

**Client:** Spectra East Inc.

**Contract:** SPEC01

**Lab code:** ACE

**Sdg no.:** Q3551

**Instrument id number:** \_\_\_\_\_

**Method:** \_\_\_\_\_

**Run number:** LB137841

**Start date:** 11/10/2025

**End date:** 11/10/2025

| Lab sample id. | Client Sample Id | d/f | Time | Parameter list                |
|----------------|------------------|-----|------|-------------------------------|
| CCB10          | CCB10            | 1   | 2247 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCV11          | CCV11            | 1   | 2330 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |
| CCB11          | CCB11            | 1   | 2334 | Ag,As,Be,Cd,Cr,Cu,Ni,Pb,Se,Zn |