

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL101425

Review By	Semsettin Yesilyurt	Review On	10/14/2025 2:02:14 PM
Supervise By	Maresh Dadoda	Supervise On	10/16/2025 4:41:37 AM
SubDirectory	VL101425	HP Acquire Method	TO15AIR
		HP Processing Method	VL101425AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2684
Initial Calibration Stds	AP2687,AP2688,AP2689
CCC	AP2687
Internal Standard/PEM	AP2684
ICV/I.BLK	AP2690
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	Sampled	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL043103.D	14 Oct 2025 08:42		SY/MD	Ok
2	VSTDICC0.5	VSTDICC0.5	VL043104.D	14 Oct 2025 09:57	not use	SY/MD	Not Ok
3	VSTDICCC010	VSTDICCC010	VL043105.D	14 Oct 2025 10:31		SY/MD	Ok,M
4	VSTDICC002	VSTDICC002	VL043106.D	14 Oct 2025 11:27		SY/MD	Ok,M
5	VSTDICC001	VSTDICC001	VL043107.D	14 Oct 2025 12:01		SY/MD	Ok,M
6	VSTDICC0.5	VSTDICC0.5	VL043108.D	14 Oct 2025 12:34		SY/MD	Ok,M
7	VSTDICC0.1	VSTDICC0.1	VL043109.D	14 Oct 2025 13:08		SY/MD	Ok,M
8	VSTDICC0.03	VSTDICC0.03	VL043110.D	14 Oct 2025 13:42		SY/MD	Ok,M
9	VSTDICC015	VSTDICC015	VL043111.D	14 Oct 2025 14:18		SY/MD	Ok,M
10	VSTDICV010	ICVVL101425	VL043112.D	14 Oct 2025 14:51		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL110625

Review By	Semsettin Yesilyurt	Review On	11/10/2025 11:08:23 AM
Supervise By	Maresh Dadoda	Supervise On	11/10/2025 1:02:12 PM
SubDirectory	VL110625	HP Acquire Method	TO15AIR
		HP Processing Method	VL101425AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2684
Initial Calibration Stds	AP2687,AP2688,AP2689
CCC	AP2687
Internal Standard/PEM	AP2684
ICV/I.BLK	AP2690
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL043154.D	06 Nov 2025 06:24		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL043155.D	06 Nov 2025 07:19		SY/MD	Ok,M
3	VL1106ABL01	VL1106ABL01	VL043156.D	06 Nov 2025 08:02		SY/MD	Ok
4	VL1106ABS01	VL1106ABS01	VL043157.D	06 Nov 2025 08:45		SY/MD	Ok,M
5	Q3555-02	151-IA-2	VL043158.D	06 Nov 2025 09:42		SY/MD	Ok,M
6	Q3555-02DUP	151-IA-2DUP	VL043159.D	06 Nov 2025 10:17		SY/MD	Ok,M
7	Q3555-01	151-IA-1	VL043160.D	06 Nov 2025 10:51		SY/MD	Ok,M
8	Q3556-01	131-IA1	VL043161.D	06 Nov 2025 11:36		SY/MD	Ok,M
9	Q3556-02	131-IA2	VL043162.D	06 Nov 2025 12:10		SY/MD	Ok,M
10	Q3557-01	141-AA-1	VL043163.D	06 Nov 2025 12:44		SY/MD	Ok,M
11	Q3557-02	141-IA-1	VL043164.D	06 Nov 2025 13:17	Need 4X	SY/MD	Dilution
12	Q3557-03	141-IA2	VL043165.D	06 Nov 2025 13:50	Need 4X	SY/MD	Dilution
13	Q3555-04DL	151-SG-2DL	VL043166.D	06 Nov 2025 14:22	4	SY/MD	Ok
14	Q3555-03DL	151-SG-1DL	VL043167.D	06 Nov 2025 14:54	3	SY/MD	Ok,M
15	Q3556-03DL	131-SG-1DL	VL043168.D	06 Nov 2025 15:27	Run @ Lower Dilution	SY/MD	Not Ok
16	Q3556-04DL	131-SG-2DL	VL043169.D	06 Nov 2025 15:59	Run @ Lower Dilution	SY/MD	Not Ok
17	Q3557-04DL	141-SG-1DL	VL043170.D	06 Nov 2025 16:31	Run @ Lower Dilution	SY/MD	Not Ok
18	Q3557-05DL	141-SG-2DL	VL043171.D	06 Nov 2025 17:03		SY/MD	Ok

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL110625

Review By	Semsettin Yesilyurt	Review On	11/10/2025 11:08:23 AM		
Supervise By	Mahesh Dadoda	Supervise On	11/10/2025 1:02:12 PM		
SubDirectory	VL110625	HP Acquire Method	TO15AIR	HP Processing Method	VL101425AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2684				
Initial Calibration Stds	AP2687,AP2688,AP2689				
CCC	AP2687				
Internal Standard/PEM	AP2684				
ICV/I.BLK	AP2690				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q3557-02DL	141-IA-1DL	VL043172.D	06 Nov 2025 17:35		SY/MD	Ok,M
20	Q3555-04	151-SG-2	VL043173.D	06 Nov 2025 18:23	Need 10X	SY/MD	Dilution
21	Q3555-03	151-SG-1	VL043174.D	06 Nov 2025 18:57	Need 10X	SY/MD	Dilution
22	Q3556-03	131-SG-1	VL043175.D	06 Nov 2025 19:30		SY/MD	Ok,M
23	Q3556-04	131-SG-2	VL043176.D	06 Nov 2025 20:03		SY/MD	Ok,M
24	Q3557-04	141-SG-1	VL043177.D	06 Nov 2025 20:35		SY/MD	Ok,M
25	Q3557-05	141-SG-2	VL043178.D	06 Nov 2025 21:08	Need 10X	SY/MD	Dilution

M : Manual Integration