

Surrogate Summary

SDG No.: Q3569

Client: Remington & Vernick Engineers

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO114870.D	PIBLK-PO114870.D	Tetrachloro-m-xyl	1	20	15.5	77		70 (60)	130 (140)
		Decachlorobiphen	1	20	16.8	84		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	16.2	81		70 (60)	130 (140)
		Decachlorobiphen	2	20	17.3	86		70 (60)	130 (140)
I.BLK-PO115014.D	PIBLK-PO115014.D	Tetrachloro-m-xyl	1	20	19.8	99		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.8	99		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	19.9	100		70 (60)	130 (140)
		Decachlorobiphen	2	20	21.7	108		70 (60)	130 (140)
PB170484BL	PB170484BL	Tetrachloro-m-xyl	1	20	23.3	117		30 (30)	150 (173)
		Decachlorobiphen	1	20	23.3	117		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	23.5	118		30 (30)	150 (173)
		Decachlorobiphen	2	20	25.5	128		30 (10)	150 (173)
PB170484BS	PB170484BS	Tetrachloro-m-xyl	1	20	24.5	123		30 (30)	150 (173)
		Decachlorobiphen	1	20	24.5	122		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	22.5	113		30 (30)	150 (173)
		Decachlorobiphen	2	20	26.8	134		30 (10)	150 (173)
PB170484BSD	PB170484BSD	Tetrachloro-m-xyl	1	20	24.3	121		30 (30)	150 (173)
		Decachlorobiphen	1	20	24.4	122		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	22.8	114		30 (30)	150 (173)
		Decachlorobiphen	2	20	26.1	130		30 (10)	150 (173)
Q3569-01	MW-2	Tetrachloro-m-xyl	1	20	24.4	122		30 (30)	150 (173)
		Decachlorobiphen	1	20	16.2	81		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	24.2	121		30 (30)	150 (173)
		Decachlorobiphen	2	20	17.4	87		30 (10)	150 (173)
Q3569-02	MW-12	Tetrachloro-m-xyl	1	20	23.9	120		30 (30)	150 (173)
		Decachlorobiphen	1	20	20.6	103		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	24.9	125		30 (30)	150 (173)
		Decachlorobiphen	2	20	22.0	110		30 (10)	150 (173)
I.BLK-PO115029.D	PIBLK-PO115029.D	Tetrachloro-m-xyl	1	20	19.5	97		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.6	98		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	19.5	97		70 (60)	130 (140)
		Decachlorobiphen	2	20	20.7	103		70 (60)	130 (140)