

Surrogate Summary

SDG No.: **Q3569**

Client: **Remington & Vernick Engineers**

Analytical Method: **8081B**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PD090647.D	PIBLK-PD090647.D	Decachlorobiphen	1	20	21.3	106		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	19.4	97		30 (41)	150 (134)
		Decachlorobiphen	2	20	20.3	102		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	20.0	100		30 (41)	150 (134)
I.BLK-PD091081.D	PIBLK-PD091081.D	Decachlorobiphen	1	20	21.9	109		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	21.8	109		30 (41)	150 (134)
		Decachlorobiphen	2	20	27.4	137		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	25.6	128		30 (41)	150 (134)
PB170485BL	PB170485BL	Decachlorobiphen	1	20	20.7	103		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	20.4	102		30 (41)	150 (134)
		Decachlorobiphen	2	20	25.6	128		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	23.6	118		30 (41)	150 (134)
PB170485BS	PB170485BS	Decachlorobiphen	1	20	21.4	107		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	20.5	103		30 (41)	150 (134)
		Decachlorobiphen	2	20	26.4	132		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	23.9	120		30 (41)	150 (134)
PB170485BSD	PB170485BSD	Decachlorobiphen	1	20	17.5	87		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	19.9	99		30 (41)	150 (134)
		Decachlorobiphen	2	20	20.7	104		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	23.1	115		30 (41)	150 (134)
Q3569-01	MW-2	Decachlorobiphen	1	20	12.7	64		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	18.0	90		30 (41)	150 (134)
		Decachlorobiphen	2	20	19.0	95		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	19.7	99		30 (41)	150 (134)
Q3569-02	MW-12	Decachlorobiphen	1	20	16.5	82		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	19.0	95		30 (41)	150 (134)
		Decachlorobiphen	2	20	81.0	405	*	30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	17.1	85		30 (41)	150 (134)
I.BLK-PD091104.D	PIBLK-PD091104.D	Decachlorobiphen	1	20	21.7	108		30 (31)	150 (149)
		Tetrachloro-m-xyl	1	20	22.0	110		30 (41)	150 (134)
		Decachlorobiphen	2	20	25.9	130		30 (31)	150 (149)
		Tetrachloro-m-xyl	2	20	24.7	124		30 (41)	150 (134)