

WORKLIST(Hardcopy Internal Chain)

WorkList Name : PB170729

WorkList ID : 193342

Department : Digestion

Date : 11-25-2025 09:56:15

Sample	Customer Sample	Matrix	Test	Preservative	Customer	Raw Sample Storage Location	Collect Date	Method
Q3716-01	MW5	Water	Metals Group4	1:1 HNO3 to pH < 2	GENV01	E22	11/24/2025	6010D
Q3716-02	MW3	Water	Metals Group4	1:1 HNO3 to pH < 2	GENV01	E22	11/24/2025	6010D
Q3716-03	MW4	Water	Metals Group4	1:1 HNO3 to pH < 2	GENV01	E22	11/24/2025	6010D
Q3716-04	MW7	Water	Metals Group4	1:1 HNO3 to pH < 2	GENV01	E22	11/24/2025	6010D
Q3718-01	Basement Kitchen (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-02	2nd draw	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-03	Adult Bathroom (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-04	2nd draw	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-05	Drinking Fountain (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-06	2nd draw	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-07	Kids Bathroom L (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-08	L 2nd draw	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-09	R (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-10	R (2nd draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-11	Main Hall (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-12	L (2nd draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-13	R (1st draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D
Q3718-14	R (2nd draw)	Water	Metals Group3	1:1 HNO3 to pH < 2	TOTA02	E61	11/24/2025	6010D

Date/Time 11/25/2025 10:20
 Raw Sample Received by: SKS met dia
 Raw Sample Relinquished by: CAP CS met

Date/Time 11/25/2025 11:20
 Raw Sample Received by: CAP CS met
 Raw Sample Relinquished by: SKS met dia