

Surrogate Summary

SDG No.: Q3743

Client: Remington & Vernick

Analytical Method: 608.3 PCB

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO114870.D	PIBLK-PO114870.D	Tetrachloro-m-xyl	1	20	15.5	77		60	140
		Decachlorobiphen	1	20	16.8	84		60	140
		Tetrachloro-m-xyl	2	20	16.2	81		60	140
		Decachlorobiphen	2	20	17.3	86		60	140
I.BLK-PO115540.D	PIBLK-PO115540.D	Tetrachloro-m-xyl	1	20	25.1	126		60	140
		Decachlorobiphen	1	20	25.5	127		60	140
		Tetrachloro-m-xyl	2	20	23.1	116		60	140
		Decachlorobiphen	2	20	24.0	120		60	140
PB170783BL	PB170783BL	Tetrachloro-m-xyl	1	20	22.4	112		60	140
		Decachlorobiphen	1	20	22.9	114		60	140
		Tetrachloro-m-xyl	2	20	22.0	110		60	140
		Decachlorobiphen	2	20	21.6	108		60	140
PB170783BS	PB170783BS	Tetrachloro-m-xyl	1	20	23.0	115		60	140
		Decachlorobiphen	1	20	23.9	119		60	140
		Tetrachloro-m-xyl	2	20	22.8	114		60	140
		Decachlorobiphen	2	20	22.2	111		60	140
PB170783BSD	PB170783BSD	Tetrachloro-m-xyl	1	20	22.7	114		60	140
		Decachlorobiphen	1	20	23.4	117		60	140
		Tetrachloro-m-xyl	2	20	22.5	113		60	140
		Decachlorobiphen	2	20	21.8	109		60	140
I.BLK-PO115550.D	PIBLK-PO115550.D	Tetrachloro-m-xyl	1	20	24.7	123		60	140
		Decachlorobiphen	1	20	25.6	128		60	140
		Tetrachloro-m-xyl	2	20	23.6	118		60	140
		Decachlorobiphen	2	20	23.7	119		60	140
I.BLK-PO115556.D	PIBLK-PO115556.D	Tetrachloro-m-xyl	1	20	24.9	125		60	140
		Decachlorobiphen	1	20	26.0	130		60	140
		Tetrachloro-m-xyl	2	20	23.6	118		60	140
		Decachlorobiphen	2	20	23.9	120		60	140
Q3743-02	TW-1125-2	Tetrachloro-m-xyl	1	20	12.1	60		60	140
		Decachlorobiphen	1	20	13.6	68		60	140
		Tetrachloro-m-xyl	2	20	12.1	61		60	140
		Decachlorobiphen	2	20	9.73	49	*	60	140
I.BLK-PO115569.D	PIBLK-PO115569.D	Tetrachloro-m-xyl	1	20	25.5	127		60	140
		Decachlorobiphen	1	20	25.7	129		60	140
		Tetrachloro-m-xyl	2	20	24.1	120		60	140
		Decachlorobiphen	2	20	23.8	119		60	140