

## Process Report

**Order ID :** q3744

| Date       | Bottle ID | Lab Sample Number | InComing Vacuum | Canister ID | Can Size | Reg ID | Reg Size | Batch Code | Certification ID                                                                                                     |
|------------|-----------|-------------------|-----------------|-------------|----------|--------|----------|------------|----------------------------------------------------------------------------------------------------------------------|
| 12/02/2025 | B2511037  | Q3744-03          | -4.7            | 10604       | 6 L      | 10550  |          | C2510005   | <b>VL043134.D</b><br>Seq :VL101725<br>TUNE : VL043130.D<br>ICAL : VL043110.D<br>CCAL : VL043131.D<br>MB : VL043132.D |
| 12/02/2025 | B2511037  | Q3744-04          | -4.3            | 10601       | 6 L      | 10511  |          | C2510005   | <b>VL043134.D</b><br>Seq :VL101725<br>TUNE : VL043130.D<br>ICAL : VL043110.D<br>CCAL : VL043131.D<br>MB : VL043132.D |
| 12/02/2025 | B2511037  | Q3744-05          | -2.7            | 10443       | 6 L      | 10648  |          | C2510005   | <b>VL043134.D</b><br>Seq :VL101725<br>TUNE : VL043130.D<br>ICAL : VL043110.D<br>CCAL : VL043131.D<br>MB : VL043132.D |
| 12/02/2025 | B2511037  | Q3744-01          | -5.1            | 10311       | 6 L      | 10185  |          | C2510005   | <b>VL043134.D</b><br>Seq :VL101725<br>TUNE : VL043130.D<br>ICAL : VL043110.D<br>CCAL : VL043131.D<br>MB : VL043132.D |
| 12/02/2025 | B2511037  | Q3744-02          | -3.5            | 10285       | 6 L      | 10707  |          | C2510005   | <b>VL043134.D</b><br>Seq :VL101725<br>TUNE : VL043130.D<br>ICAL : VL043110.D<br>CCAL : VL043131.D<br>MB : VL043132.D |