

Metals

- 7 -

LABORATORY CONTROL SAMPLE SUMMARY

Client: Tetra Tech, EMI

SDG No.: Q3764

Contract: TETR16

Lab Code: ACE

Analyte	Units	True Value	Result	C	% Recovery	Acceptance Limits	M
PB170827BS							
Aluminum	ug/L	1000	978		98	80 - 120	P
Antimony	ug/L	400	413		103	80 - 120	P
Arsenic	ug/L	400	408		102	80 - 120	P
Barium	ug/L	100	92.7		93	80 - 120	P
Beryllium	ug/L	100	97.2		97	80 - 120	P
Cadmium	ug/L	100	98.7		99	80 - 120	P
Calcium	ug/L	500	495	J	99	80 - 120	P
Chromium	ug/L	200	202		101	80 - 120	P
Cobalt	ug/L	100	98.6		99	80 - 120	P
Copper	ug/L	150	158		105	80 - 120	P
Iron	ug/L	1500	1460		97	80 - 120	P
Lead	ug/L	500	486		97	80 - 120	P
Magnesium	ug/L	1000	983	J	98	80 - 120	P
Manganese	ug/L	100	95.6		96	80 - 120	P
Nickel	ug/L	250	248		99	80 - 120	P
Potassium	ug/L	5000	4910		98	80 - 120	P
Selenium	ug/L	1000	1060		106	80 - 120	P
Silver	ug/L	37.5	38.3		102	80 - 120	P
Sodium	ug/L	1500	1500		100	80 - 120	P
Thallium	ug/L	1000	1010		101	80 - 120	P
Vanadium	ug/L	150	146		97	80 - 120	P
Zinc	ug/L	100	105		105	80 - 120	P