

Metals

- 7 -

LABORATORY CONTROL SAMPLE SUMMARY

Client: Tetra Tech, EMI

SDG No.: Q3766

Contract: TETR16

Lab Code: ACE

Analyte	Units	True Value	Result	C	% Recovery	Acceptance Limits	M
PB170842BS							
Aluminum	ug/L	1000	920		92	80 - 120	P
Antimony	ug/L	400	392		98	80 - 120	P
Arsenic	ug/L	400	380		95	80 - 120	P
Barium	ug/L	100	89.1		89	80 - 120	P
Beryllium	ug/L	100	91.4		91	80 - 120	P
Cadmium	ug/L	100	93.1		93	80 - 120	P
Calcium	ug/L	500	470	J	94	80 - 120	P
Chromium	ug/L	200	191		96	80 - 120	P
Cobalt	ug/L	100	93.7		94	80 - 120	P
Copper	ug/L	150	149		99	80 - 120	P
Iron	ug/L	1500	1390		93	80 - 120	P
Lead	ug/L	500	460		92	80 - 120	P
Magnesium	ug/L	1000	934	J	93	80 - 120	P
Manganese	ug/L	100	92.2		92	80 - 120	P
Nickel	ug/L	250	235		94	80 - 120	P
Potassium	ug/L	5000	4530		91	80 - 120	P
Selenium	ug/L	1000	982		98	80 - 120	P
Silver	ug/L	37.5	35.7		95	80 - 120	P
Sodium	ug/L	1500	1380		92	80 - 120	P
Thallium	ug/L	1000	923		92	80 - 120	P
Vanadium	ug/L	150	139		93	80 - 120	P
Zinc	ug/L	100	100		100	80 - 120	P