

Surrogate Summary

SW-846

SDG No.: Q3834

Client: Woodard & Curran

Analytical Method: 8270E

Lab Sample ID	Client ID	Parameter	Spike (PPM)	Result (PPM)	Recovery (%)	Qual	Limits (%)	
							Low	High
PB170902BL	PB170902BL	2-Fluorophenol	150	130	87		30 (10)	130 (105)
		Phenol-d6	150	121	81		30 (10)	130 (103)
		Nitrobenzene-d5	100	87.0	87		30 (10)	130 (108)
		2-Fluorobiphenyl	100	89.3	89		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	147	98		30 (10)	130 (116)
		Terphenyl-d14	100	91.2	91		30 (10)	130 (109)
PB170902BS	PB170902BS	2-Fluorophenol	150	117	78		30 (10)	130 (105)
		Phenol-d6	150	113	75		30 (10)	130 (103)
		Nitrobenzene-d5	100	76.2	76		30 (10)	130 (108)
		2-Fluorobiphenyl	100	79.0	79		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	134	89		30 (10)	130 (116)
		Terphenyl-d14	100	84.2	84		30 (10)	130 (109)
Q3831-01MS	PL-01-12102025MS	2-Fluorophenol	150	67.7	45		30 (10)	130 (105)
		Phenol-d6	150	66.2	44		30 (10)	130 (103)
		Nitrobenzene-d5	100	45.4	45		30 (10)	130 (108)
		2-Fluorobiphenyl	100	45.2	45		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	74.3	50		30 (10)	130 (116)
		Terphenyl-d14	100	46.4	46		30 (10)	130 (109)
Q3831-01MSD	PL-01-12102025MSD	2-Fluorophenol	150	64.0	43		30 (10)	130 (105)
		Phenol-d6	150	62.2	41		30 (10)	130 (103)
		Nitrobenzene-d5	100	43.2	43		30 (10)	130 (108)
		2-Fluorobiphenyl	100	43.3	43		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	71.8	48		30 (10)	130 (116)
		Terphenyl-d14	100	45.7	46		30 (10)	130 (109)
Q3834-29	B127(5.5-8)120925	2-Fluorophenol	150	63.2	42		30 (10)	130 (105)
		Phenol-d6	150	58.3	39		30 (10)	130 (103)
		Nitrobenzene-d5	100	39.7	40		30 (10)	130 (108)
		2-Fluorobiphenyl	100	36.4	36		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	63.9	43		30 (10)	130 (116)
		Terphenyl-d14	100	35.6	36		30 (10)	130 (109)
Q3834-30	B127(8-10)120925	2-Fluorophenol	150	71.3	48		30 (10)	130 (105)
		Phenol-d6	150	67.3	45		30 (10)	130 (103)
		Nitrobenzene-d5	100	46.5	47		30 (10)	130 (108)
		2-Fluorobiphenyl	100	45.5	46		30 (10)	130 (108)
		2,4,6-Tribromophenol	150	76.7	51		30 (10)	130 (116)
		Terphenyl-d14	100	42.9	43		30 (10)	130 (109)