

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_D\Data\PD051625\
 Data File : PD088578.D
 Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
 Acq On : 16 May 2025 10:40
 Operator : AR\AJ
 Sample : Q1994-12
 Misc :
 ALS Vial : 5 Sample Multiplier: 1

Instrument :
 ECD_D
 ClientSampleId :
 PEST-GPC2-BLANK-SPIKE

Manual Integrations
 APPROVED

Reviewed By :Abdul Mirza 05/22/2025
 Supervised By :mohammad ahmed 05/23/2025

Integration File signal 1: autoint1.e
 Integration File signal 2: autoint2.e
 Quant Time: May 22 14:20:07 2025
 Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_D\Method\PD050925CLP.M
 Quant Title : GC Extractables
 QLast Update : Mon May 12 12:22:29 2025
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 1 µl
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
 Signal #1 Info : 30M x 0.32mm x0.5 Signal #2 Info : 30M x 0.32mm x 0.25µm

Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml
----------	------	------	--------	--------	-------	-------

System Monitoring Compounds

Target Compounds

3) MA	gamma-BHC...	4.338	3.752	96973007	477.0E6	21.092	20.672m
4) MA	Heptachlor	4.937	4.106	92406401	450.6E6	19.535	18.900
5) MB	Aldrin	5.279	4.392	97059176	478.4E6	21.757	21.302
13) MA	Dieldrin	6.354	5.535	175.8E6	862.7E6	41.119	40.630m
14) MA	Endrin	6.582	5.811	139.9E6	734.8E6	39.031	38.472m
17) MA	4,4'-DDT	7.029	6.205	105.3E6	607.6E6	33.627	34.456m

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_D\Data\PD051625\
Data File : PD088578.D
Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
Acq On : 16 May 2025 10:40
Operator : AR\AJ
Sample : Q1994-12
Misc :
ALS Vial : 5 Sample Multiplier: 1

Integration File signal 1: autoint1.e
Integration File signal 2: autoint2.e
Quant Time: May 22 14:20:07 2025
Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_D\Method\PD050925CLP.M
Quant Title : GC Extractables
QLast Update : Mon May 12 12:22:29 2025
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 1 µl
Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
Signal #1 Info : 30M x 0.32mm x 0.5 Signal #2 Info : 30M x 0.32mm x 0.25µm

