

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD\_D\Data\PD072924\  
 Data File : PD084116.D  
 Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch  
 Acq On : 29 Jul 2024 11:07  
 Operator : AR\AJ  
 Sample : PB162310BL  
 Misc :  
 ALS Vial : 9 Sample Multiplier: 1

**Instrument :**  
 ECD\_D  
**ClientSampleId :**  
 PB162310BL

**Manual Integrations**  
**APPROVED**

Reviewed By :Abdul Mirza 07/30/2024  
 Supervised By :Ankita Jodhani 07/30/2024

Integration File signal 1: autoint1.e  
 Integration File signal 2: autoint2.e  
 Quant Time: Jul 30 02:37:34 2024  
 Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD\_D\Method\PD071524.M  
 Quant Title : GC Extractables  
 QLast Update : Tue Jul 16 04:22:06 2024  
 Response via : Initial Calibration  
 Integrator: ChemStation

Volume Inj. : 1 µl  
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2  
 Signal #1 Info : 30M x 0.32mm x0.5 Signal #2 Info : 30M x 0.32mm x 0.25µm

Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml
-----						
System Monitoring Compounds						
1) SA Tetrachlo...	3.567	2.768	17083467	91204178	16.708m	15.863
28) SA Decachlor...	9.127	7.896	37967947	117.8E6	16.632	18.351

Target Compounds

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD\_D\Data\PD072924\  
 Data File : PD084116.D  
 Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch  
 Acq On : 29 Jul 2024 11:07  
 Operator : AR\AJ  
 Sample : PB162310BL  
 Misc :  
 ALS Vial : 9 Sample Multiplier: 1

**Instrument :**

ECD\_D

**ClientSampleId :**

PB162310BL

**Manual Integrations**

**APPROVED**

Reviewed By :Abdul Mirza 07/30/2024

Supervised By :Ankita Jodhani 07/30/2024

Integration File signal 1: autoint1.e  
 Integration File signal 2: autoint2.e  
 Quant Time: Jul 30 02:37:34 2024  
 Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD\_D\Method\PD071524.M  
 Quant Title : GC Extractables  
 QLast Update : Tue Jul 16 04:22:06 2024  
 Response via : Initial Calibration  
 Integrator: ChemStation

Volume Inj. : 1 µl  
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2  
 Signal #1 Info : 30M x 0.32mm x0.5 Signal #2 Info : 30M x 0.32mm x 0.25µm

