

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_D\Data\PD122319\
 Data File : PD056745.D
 Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
 Acq On : 23 Dec 2019 18:04
 Operator : SG\AJ
 Sample : K6322-13MS
 Misc :
 ALS Vial : 32 Sample Multiplier: 1

Instrument :
 ECD_D
 ClientSampleId :

Integration File signal 1: autoint1.e
 Integration File signal 2: autoint2.e
 Quant Time: Dec 24 00:40:06 2019
 Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_D\Method\PD120219CLP.M
 Quant Title : GC Extractables
 QLast Update : Tue Dec 03 08:51:13 2019
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 1 µl
 Signal #1 Phase : ZB-MR2 Signal #2 Phase: ZB-MR1
 Signal #1 Info : 30M x 0.32mm x 0.2 Signal #2 Info : 30M x 0.32mm x 0.50µm

	Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml

System Monitoring Compounds							
1)	SA Tetrachlo...	3.286	3.961	10013107	17931687	12.090	14.915
27)	SA Decachlor...	7.958	8.996	13901584	20262680	15.395	16.841m
Target Compounds							
3)	MA gamma-BHC...	3.988	4.630	49740300	86919931	44.295	53.859
4)	MA Heptachlor	4.273	5.136	48144797	85909259	47.767	52.712
5)	MB Aldrin	4.512	5.439	46648988	85596690	42.955	52.092
13)	MA Dieldrin	5.519	6.435	103.0E6	160.0E6	101.395	101.830
14)	MA Endrin	5.771	6.651	79831141	138.0E6	95.527m	109.019
17)	MA 4,4'-DDT	6.142	7.071	77394749	124.9E6	109.992	105.646m

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_D\Data\PD122319\
Data File : PD056745.D
Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
Acq On : 23 Dec 2019 18:04
Operator : SG\AJ
Sample : K6322-13MS
Misc :
ALS Vial : 32 Sample Multiplier: 1

Instrument :
ECD_D
ClientSampleId :

Integration File signal 1: autoint1.e
Integration File signal 2: autoint2.e
Quant Time: Dec 24 00:40:06 2019
Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_D\Method\PD120219CLP.M
Quant Title : GC Extractables
QLast Update : Tue Dec 03 08:51:13 2019
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 1 µl
Signal #1 Phase : ZB-MR2 Signal #2 Phase: ZB-MR1
Signal #1 Info : 30M x 0.32mm x0.2 Signal #2 Info : 30M x 0.32mm x 0.50µm

