

Data Path : P:\HPCHEM1\ECD_0\Data\P0021918\
 Data File : P0042471.D
 Signal(s) : Signal #1: ECD1A.CH Signal #2: ECD2B.CH
 Acq On : 19 Feb 2018 11:15
 Operator : SM/UA
 Sample : AR1254ICC500
 Misc :
 ALS Vial : 12 Sample Multiplier: 1

Instrument :
 ECD_O
 ClientSampleId :
 AR1254ICC500

Integration File signal 1: autoint1.e
 Integration File signal 2: autoint2.e
 Quant Time: Feb 19 12:50:19 2018
 Quant Method : P:\HPCHEM1\ECD_0\P0021918.M
 Quant Title : GC EXTRACTABLES
 QLast Update : Mon Feb 19 12:39:21 2018
 Response via : Initial Calibration
 Integrator: ChemStation 6890 Scale Mode: Small noise peaks clipped

Volume Inj. : 2 µl
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
 Signal #1 Info : 30Mx0.32mmx 0.50µ Signal #2 Info : 30M x 0.32mm x 0.25µm

Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml

System Monitoring Compounds						
1) SA Tetrachlo...	4.093	3.322	42037339	53362306	50.000	50.000
2) SA Decachlor...	9.500	8.144	31845252	39047201	50.000	50.000
Target Compounds						
26) L6 AR-1254-1	6.278	5.145	21031595	26426879	500.000	500.000
27) L6 AR-1254-2	6.551	5.291	12851691	22594925	500.000	500.000
28) L6 AR-1254-3	6.635	5.680	23336395	36470422	500.000	500.000
29) L6 AR-1254-4	6.915	5.905	19007305	25047050	500.000	500.000
30) L6 AR-1254-5	7.180	6.209	13670704	16401817	500.000	500.000

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : P:\HPCHEM1\ECD_0\Data\P0021918\
Data File : P0042471.D
Signal(s) : Signal #1: ECD1A.CH Signal #2: ECD2B.CH
Acq On : 19 Feb 2018 11:15
Operator : SM/UA
Sample : AR1254ICC500
Misc :
ALS Vial : 12 Sample Multiplier: 1

Instrument :
ECD_O
ClientSampleId :
AR1254ICC500

Integration File signal 1: autoint1.e
Integration File signal 2: autoint2.e
Quant Time: Feb 19 12:50:19 2018
Quant Method : P:\HPCHEM1\ECD_0\P0021918.M
Quant Title : GC EXTRACTABLES
QLast Update : Mon Feb 19 12:39:21 2018
Response via : Initial Calibration
Integrator: ChemStation 6890 Scale Mode: Small noise peaks clipped

Volume Inj. : 2 µl
Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
Signal #1 Info : 30Mx0.32mmx 0.50µ Signal #2 Info : 30M x 0.32mm x 0.25µm

