

Data Path : P:\HPCHEM1\ECD_0\Data\P0021918\
 Data File : P0042472.D
 Signal(s) : Signal #1: ECD1A.CH Signal #2: ECD2B.CH
 Acq On : 19 Feb 2018 11:30
 Operator : SM/UA
 Sample : AR1262ICC500
 Misc :
 ALS Vial : 13 Sample Multiplier: 1

Instrument :
 ECD_O
 ClientSampleId :
 AR1262ICC500

Integration File signal 1: autoint1.e
 Integration File signal 2: autoint2.e
 Quant Time: Feb 19 13:00:37 2018
 Quant Method : P:\HPCHEM1\ECD_0\P0021918.M
 Quant Title : GC EXTRACTABLES
 QLast Update : Mon Feb 19 12:49:13 2018
 Response via : Initial Calibration
 Integrator: ChemStation 6890 Scale Mode: Small noise peaks clipped

Volume Inj. : 2 µl
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
 Signal #1 Info : 30Mx0.32mmx 0.50µ Signal #2 Info : 30M x 0.32mm x 0.25µm

Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml

System Monitoring Compounds						
1) SA Tetrachlo...	4.086	3.322	43113384	56066807	50.000	50.000
2) SA Decachlor...	9.491	8.140	39311783	49824074	50.000	50.000
Target Compounds						
36) L8 AR-1262-1	6.791	5.991	13497016	28535683	500.000	500.000
37) L8 AR-1262-2	7.044	6.137	18013818	24960594	500.000	500.000
38) L8 AR-1262-3	7.317	6.345	16582310	46260560	500.000	500.000
39) L8 AR-1262-4	7.618	6.596	22978960	34471080	500.000	500.000
40) L8 AR-1262-5	7.918	6.839	52822307	82592922	500.000	500.000

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : P:\HPCHEM1\ECD_0\Data\P0021918\
Data File : P0042472.D
Signal(s) : Signal #1: ECD1A.CH Signal #2: ECD2B.CH
Acq On : 19 Feb 2018 11:30
Operator : SM/UA
Sample : AR1262ICC500
Misc :
ALS Vial : 13 Sample Multiplier: 1

Instrument :
ECD_O
ClientSampleId :
AR1262ICC500

Integration File signal 1: autoint1.e
Integration File signal 2: autoint2.e
Quant Time: Feb 19 13:00:37 2018
Quant Method : P:\HPCHEM1\ECD_0\P0021918.M
Quant Title : GC EXTRACTABLES
QLast Update : Mon Feb 19 12:49:13 2018
Response via : Initial Calibration
Integrator: ChemStation 6890 Scale Mode: Small noise peaks clipped

Volume Inj. : 2 µl
Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
Signal #1 Info : 30Mx0.32mmx 0.50µ Signal #2 Info : 30M x 0.32mm x 0.25µm

