

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_R\Data\PR022423\
 Data File : PR059773.D
 Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
 Acq On : 25 Feb 2023 03:07
 Operator : YP\AJ
 Sample : 01627-08
 Misc : AR1232 LOQ 50 PPB
 ALS Vial : 50 Sample Multiplier: 1

Instrument :
 ECD_R
 ClientSampleId :
 LOQ-WATER-02-QT1-2023

Integration File signal 1: autoint1.e
 Integration File signal 2: autoint2.e
 Quant Time: Feb 26 04:55:13 2023
 Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_R\Method\PR022423.M
 Quant Title : GC EXTRACTABLES
 QLast Update : Fri Feb 24 12:41:32 2023
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 2 µl
 Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
 Signal #1 Info : 30Mx0.32mmx 0.50µ Signal #2 Info : 30M x 0.32mm x 0.25µm

Compound	RT#1	RT#2	Resp#1	Resp#2	ng/ml	ng/ml

System Monitoring Compounds						
1) SA Tetrachlo...	3.690	2.906	56969308	82366256	19.068	19.029
2) SA Decachlor...	9.658	8.137	52198641	76006799	22.757	21.545
Target Compounds						
11) L3 AR-1232-1	4.081	3.291	3780694	4549616	54.357	49.176
12) L3 AR-1232-2	4.636	4.034	1716166	4823225	52.753	54.722
13) L3 AR-1232-3	4.947	4.212	3195106	1813830	54.591	40.600 #
14) L3 AR-1232-4	5.114	4.303	1329343	2062256	44.871	51.347
15) L3 AR-1232-5	5.219	4.478	1307812	2273426	58.199	49.742

(f)=RT Delta > 1/2 Window (#)=Amounts differ by > 25% (m)=manual int.

Data Path : Z:\pestpcbsrv\HPCHEM1\ECD_R\Data\PR022423\
Data File : PR059773.D
Signal(s) : Signal #1: ECD1A.ch Signal #2: ECD2B.ch
Acq On : 25 Feb 2023 03:07
Operator : YP\AJ
Sample : 01627-08
Misc : AR1232 LOQ 50 PPB
ALS Vial : 50 Sample Multiplier: 1

Instrument :
ECD_R
ClientSampleId :
LOQ-WATER-02-QT1-2023

Integration File signal 1: autoint1.e
Integration File signal 2: autoint2.e
Quant Time: Feb 26 04:55:13 2023
Quant Method : Z:\pestpcbsrv\HPCHEM1\ECD_R\Method\PR022423.M
Quant Title : GC EXTRACTABLES
QLast Update : Fri Feb 24 12:41:32 2023
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 2 µl
Signal #1 Phase : ZB-MR1 Signal #2 Phase: ZB-MR2
Signal #1 Info : 30Mx0.32mmx 0.50µm Signal #2 Info : 30M x 0.32mm x 0.25µm

