



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: jignesh
Date: 12/27/2018

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 12/26/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 12/27/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1

QC:LB100044

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
J6195-11	OR-01-122618-A	1	1.14	9.96	8.96	88.7
J6195-13	OR-01-122618-B	2	1.15	9.91	9.11	90.9
J6533-23	09017	3	1.00	2.00	2.00	100.0
J6550-01	P001-SS221-0206-01	4	1.16	9.51	8.56	88.6
J6550-02	P001-SS222-0002-01	5	1.14	9.58	8.4	86.0
J6550-03	P001-SS223-0002-01	6	1.16	9.7	7.99	80.0
J6550-04	P001-SS226-0002-01	7	1.18	9.74	8.55	86.1
J6550-05	P001-SS226-0206-01	8	1.15	9.96	9.24	91.8
J6550-06	P001-SS227-0002-01	9	1.13	9.58	8.28	84.6
J6550-07	J6550-06MS	10	1.13	9.58	8.28	84.6
J6550-08	J6550-06MSD	11	1.13	9.58	8.28	84.6
J6550-09	P001-SS227-0002-02	12	1.17	9.55	8.22	84.1
J6550-10	P001-SS227-0206-01	13	1.14	9.69	8.8	89.6
J6550-11	P001-SS227-0612-01	14	1.13	9.55	8.25	84.6
J6550-12	P002-SS007-0002-01	15	1.11	9.74	8.18	81.9
J6550-13	P002-SS008-0206-01	16	1.11	9.59	8.68	89.3
J6550-14	P002-SS010-0002-01	17	1.13	9.98	8.81	86.8
J6550-15	P002-SS011-0002-01	18	1.12	9.97	8.48	83.2
J6550-16	P002-SS012-0002-01	19	1.14	9.68	8.03	80.7
J6552-01	P001-SS204-0002-01	20	1.12	9.73	8.29	83.3
J6552-02	P001-SS204-0206-01	21	1.13	9.66	8.85	90.5
J6552-03	P001-SS204-0612-01	22	1.13	9.97	9.22	91.5
J6552-04	P001-SS205-0002-01	23	1.11	9.88	8.97	89.6
J6552-05	P001-SS206-0206-01	24	1.14	9.88	9.03	90.3
J6552-06	P001-SS207-0206-01	25	1.12	9.98	9.15	90.6
J6552-07	P001-SS208-0002-01	26	1.12	9.95	7.83	76.0
J6552-08	P001-SS209-0002-01	27	1.1	9.6	8.57	87.9
J6552-09	P001-SS210-0002-01	28	1.12	9.74	8.45	85.0
J6552-10	P001-SS211-0002-01	29	1.13	9.83	8.1	80.1
J6552-11	P001-SS213-0002-01	30	1.13	9.97	8.45	82.8
J6552-12	P001-SS214-0206-01	31	1.17	9.76	7.89	78.2
J6552-13	P001-SS215-0002-01	32	1.16	9.78	9.18	93.0
J6552-14	P001-SS216-0206-01	33	1.12	9.71	8.87	90.2
J6552-15	P001-SS217-0206-01	34	1.11	9.96	9.29	92.4
J6552-16	P001-SS218-0002-01	35	1.14	9.51	8.17	84.0
J6552-17	P001-SS219-0002-01	36	1.13	9.65	7.81	78.4
J6552-18	P001-SS220-0002-01	37	1.13	9.55	8.19	83.8
J6552-19	P001-SS220-1824-01	38	1.13	9.71	8.9	90.6
J6553-01	P001-SS184-0206-01	39	1.14	9.56	8.61	88.7



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: jignesh
Date: 12/27/2018

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 12/26/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 12/27/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1

QC:LB100044

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
J6553-02	J6553-01MS	40	1.14	9.56	8.61	88.7
J6553-03	J6553-01MSD	41	1.14	9.56	8.61	88.7
J6553-04	P001-SS185-0002-01	42	1.16	9.96	8.48	83.2
J6553-05	P001-SS187-0206-01	43	1.1	9.81	8.5	85.0
J6553-06	P001-SS187-1218-01	44	1.11	9.6	8.94	92.2
J6553-07	P001-SS188-0206-01	45	1.11	9.88	9.02	90.2
J6553-08	P001-SS189-0206-01	46	1.12	9.78	8.96	90.5
J6553-09	P001-SS190-0002-01	47	1.12	9.98	8.37	81.8
J6553-10	P001-SS190-0612-01	48	1.1	9.98	9.22	91.4
J6553-11	P001-SS191-0002-01	49	1.13	9.74	8.43	84.8
J6553-12	P001-SS192-1218-01	50	1.1	9.71	8.83	89.8
J6553-13	P001-SS193-0612-01	51	1.12	9.57	8.95	92.7
J6553-14	J6553-13MS	52	1.12	9.57	8.95	92.7
J6553-15	J6553-13MSD	53	1.12	9.57	8.95	92.7
J6553-16	P001-SS196-1218-01	54	1.13	9.51	8.88	92.5
J6553-17	P001-SS199-0206-01	55	1.13	9.93	9.22	91.9
J6553-18	P001-SS200-0206-01	56	1.11	9.93	9.2	91.7
J6553-19	P001-SS201-0206-01	57	1.13	9.52	8.67	89.9
J6553-20	P001-SS202-0002-01	58	1.13	9.97	8.81	86.9
J6553-21	P001-SS203-0002-01	59	1.13	9.63	7.78	78.2
J6553-22	P001-SS203-1218-01	60	1.13	9.96	9.06	89.8
J6553-23	P001-SS203-1824-01	61	1.12	9.89	9.15	91.6
J6554-01	P001-SS168-0206-01	62	1.12	9.64	8.85	90.7
J6554-02	P001-SS169-0002-01	63	1.12	9.65	7.85	78.9
J6554-03	P001-SS169-0612-01	64	1.12	9.68	8.26	83.4
J6554-04	P001-SS170-0206-01	65	1.11	9.72	9.00	91.6
J6554-05	P001-SS171-0612-01	66	1.13	9.83	9.18	92.5
J6554-06	P001-SS172-0002-01	67	1.11	9.58	8.49	87.1
J6554-07	P001-SS173-0002-01	68	1.12	9.69	8.38	84.7
J6554-08	P001-SS173-0206-01	69	1.1	9.96	9.26	92.1
J6554-09	J6554-08MS	70	1.1	9.96	9.26	92.1
J6554-10	J6554-08MSD	71	1.1	9.96	9.26	92.1
J6554-11	P001-SS174-0206-01	72	1.12	9.96	9.31	92.6
J6554-12	P001-SS175-0206-01	73	1.11	9.91	8.48	83.7
J6554-13	P001-SS176-0002-01	74	1.15	9.61	7.84	79.1
J6554-14	P001-SS177-0206-01	75	1.14	9.83	8.71	87.1
J6554-15	P001-SS178-0002-01	76	1.14	9.98	8.55	83.8
J6554-16	P001-SS179-0002-01	77	1.12	9.7	8.58	86.9
J6554-17	P001-SS179-0612-01	78	1.12	9.54	8.81	91.3



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: jignesh
Date: 12/27/2018

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 12/26/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 12/27/2018
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1

QC:LB100044

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
J6554-18	P001-SS180-0002-01	79	1.14	9.67	8.74	89.1
J6554-19	P001-SS181-0002-01	80	1.13	9.69	8.76	89.1
J6554-20	P001-SS183-0206-01	81	1.13	9.56	8.86	91.7
J6554-21	P001-SS183-1218-01	82	1.15	9.87	9.12	91.4
J6554-22	J6554-21MS	83	1.15	9.87	9.12	91.4
J6554-23	J6554-21MSD	84	1.15	9.87	9.12	91.4
J6568-01	SS-04	85	1.15	9.97	7.57	72.8
J6568-02	SS-03	86	1.13	9.54	7.8	79.3
J6568-03	SS-01	87	1.15	9.78	8.15	81.1
J6568-04	SS-05	88	1.11	9.52	7.65	77.8
J6568-05	SS-02	89	1.18	9.68	7.69	76.6
J6571-01	50072	90	1.00	2.00	2.00	100.0
J6571-02	10827	91	1.00	2.00	2.00	100.0
J6578-01	CF-122618-A	92	1.00	2.00	2.00	100.0
J6578-02	CF-122618-B	93	1.00	2.00	2.00	100.0
J6579-01	TN-122618-C	94	1.17	9.92	8.62	85.1
J6579-02	TN-122618-C	95	1.17	9.92	8.62	85.1
J6580-01	344	96	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$