



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/4/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:35
In Date: 02/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 02/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: % SOLIDS OVEN

QC:LB100613

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1322-01	MGAF52	1	1.18	9.96	7.74	74.7
K1322-02	MGAF53	2	1.13	9.53	7.5	75.8
K1322-03	MGAF54	3	1.18	9.68	7.84	78.4
K1322-04	MGAF55	4	1.12	9.53	7.69	78.1
K1322-05	MGAF56	5	1.18	9.96	8.08	78.6
K1322-06	MGAF57	6	1.16	9.71	7.93	79.2
K1322-07	MGAF58	7	1.14	9.91	8.06	78.9
K1322-08	MGAF59	8	1.18	9.64	7.98	80.4
K1322-09	MGAF59D	9	1.18	9.64	7.98	80.4
K1322-10	MGAF59S	10	1.18	9.64	7.98	80.4
K1322-11	MGAF60	11	1.18	9.98	8.94	88.2
K1322-12	MGAF61	12	1.1	9.75	8.41	84.5
K1322-13	MGAF62	13	1.14	9.62	7.86	79.2
K1322-14	MGAF63	14	1.16	9.5	7.67	78.1
K1322-15	MGAF64	15	1.17	9.97	8.7	85.6
K1322-16	MGAF65	16	1.13	9.89	8.2	80.7
K1322-17	MGAF66	17	1.19	9.97	8.59	84.3
K1322-18	MGAF67	18	1.14	9.96	8.43	82.7
K1322-19	MGAF68	19	1.13	9.85	8.29	82.1
K1322-20	MGAF69	20	1.16	9.82	7.92	78.1
K1322-21	MGAF70	21	1.13	9.86	8.82	88.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$