



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:15
In Date: 02/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:48
Out Date: 02/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: % SOLIDS OVEN

QC:LB100627

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1336-01	MGAFS8	1	1.13	9.6	7.34	73.3
K1336-02	MGAFS9	2	1.18	9.97	8.43	82.5
K1336-03	MGAF0	3	1.16	9.55	7.38	74.1
K1336-04	MGAF0D	4	1.16	9.55	7.38	74.1
K1336-05	MGAF0S	5	1.16	9.55	7.38	74.1
K1336-06	MGAF1	6	1.13	9.57	8.47	87.0
K1336-07	MGAF2	7	1.14	9.76	8.08	80.5
K1336-08	MGAF3	8	1.12	9.64	8.67	88.6
K1336-09	MGAF4	9	1.17	9.67	8.5	86.2
K1336-10	MGAF5	10	1.18	9.96	8.57	84.2
K1336-11	MGAF6	11	1.19	9.95	8.51	83.6
K1336-12	MGAF7	12	1.19	9.72	8.18	81.9
K1336-13	MGAF8	13	1.11	9.96	7.95	77.3
K1336-14	MGAF9	14	1.17	9.81	6.74	64.5
K1336-15	MGAFW0	15	1.13	9.53	7.46	75.4
K1336-16	MGAFW1	16	1.19	9.97	8.1	78.7
K1336-17	MGAFW2	17	1.14	9.71	7.49	74.1
K1336-18	MGAFW3	18	1.17	9.98	7.47	71.5
K1336-19	MGAFW4	19	1.16	9.5	7.08	71.0
K1336-20	MGAFW5	20	1.13	9.64	7.57	75.7
K1336-21	MGAFW6	21	1.11	9.91	8.24	81.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$