



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/6/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:10
In Date: 02/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100656

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1342-01	MGAG05	1	1.17	9.98	8.57	84.0
K1342-02	MGAG06	2	1.12	9.93	8.02	78.3
K1342-03	MGAG07	3	1.19	9.98	8.95	88.3
K1342-04	MGAG08	4	1.13	9.51	8.53	88.3
K1342-05	MGAG09	5	1.14	9.72	8.09	81.0
K1342-06	MGAG10	6	1.13	9.66	8.58	87.3
K1342-07	MGAG10D	7	1.13	9.66	8.58	87.3
K1342-08	MGAG10S	8	1.13	9.66	8.58	87.3
K1342-09	MGAG11	9	1.18	9.8	8.16	81.0
K1342-10	MGAG12	10	1.13	9.85	8.05	79.4
K1342-11	MGAG13	11	1.19	9.6	7.57	75.9
K1342-12	MGAG14	12	1.16	9.98	8.47	82.9
K1342-13	MGAG15	13	1.1	9.98	9.17	90.9
K1342-14	MGAG16	14	1.19	9.82	7.79	76.5
K1342-15	MGAG17	15	1.17	9.82	7.78	76.4
K1342-16	MGAG18	16	1.17	9.65	9.11	93.6
K1342-17	MGAG19	17	1.12	9.71	9.16	93.6
K1342-18	MGAG20	18	1.18	9.96	6.73	63.2
K1342-19	MGAG21	19	1.13	9.83	7.97	78.6
K1342-20	MGAG22	20	1.15	9.94	8.05	78.5
K1342-21	MGAG23	21	1.14	9.97	8.37	81.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$