



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/6/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:00
In Date: 02/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:41
Out Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100658

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1351-01	MGAG43	1	1.15	9.95	7.46	71.7
K1351-02	MGAG44	2	1.16	9.9	8.00	78.3
K1351-03	MGAG45	3	1.19	9.93	7.77	75.3
K1351-04	MGAG46	4	1.16	9.94	7.59	73.2
K1351-05	MGAG47	5	1.11	9.82	8.00	79.1
K1351-06	MGAG48	6	1.14	9.96	8.25	80.6
K1351-07	MGAG49	7	1.14	9.78	7.98	79.2
K1351-08	MGAG50	8	1.18	9.98	8.88	87.5
K1351-09	MGAG50D	9	1.18	9.98	8.88	87.5
K1351-10	MGAG50S	10	1.18	9.98	8.88	87.5
K1351-11	MGAG51	11	1.16	9.6	8.08	82.0
K1351-12	MGAG52	12	1.1	9.77	8.2	81.9
K1351-13	MGAG53	13	1.19	9.77	7.92	78.4
K1351-14	MGAG54	14	1.17	9.92	8.43	83.0
K1351-15	MGAG55	15	1.13	9.67	7.17	70.7
K1351-16	MGAG56	16	1.16	9.97	8.58	84.2
K1351-17	MGAG57	17	1.16	9.63	7.72	77.4
K1351-18	MGAG58	18	1.16	9.96	8.53	83.7
K1351-19	MGAG59	19	1.17	9.93	7.94	77.3
K1351-20	MGAG60	20	1.11	9.86	7.21	69.7
K1351-21	MGAG61	21	1.13	9.9	7.43	71.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$