



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/6/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:25
In Date: 02/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100662

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1357-01	MGAGA0	1	1.12	9.83	8.34	82.9
K1357-02	MGAGA1	2	1.18	9.96	8.38	82.0
K1357-03	MGAGA2	3	1.14	9.91	7.24	69.6
K1357-04	MGAGA3	4	1.11	9.96	8.63	85.0
K1357-05	MGAGA4	5	1.18	9.97	8.01	77.7
K1357-06	MGAGA5	6	1.18	9.95	8.21	80.2
K1357-07	MGAGA6	7	1.16	9.71	8.12	81.4
K1357-08	MGAGA7	8	1.12	9.79	8.32	83.0
K1357-09	MGAGA8	9	1.14	9.98	9.02	89.1
K1357-10	MGAGA9	10	1.16	9.94	6.74	63.6
K1357-11	MGAGB0	11	1.17	9.8	7.78	76.6
K1357-12	MGAGB0D	12	1.17	9.8	7.78	76.6
K1357-13	MGAGB0S	13	1.17	9.8	7.78	76.6
K1357-14	MGAGB1	14	1.19	9.98	8.08	78.4
K1357-15	MGAGB2	15	1.1	9.97	7.69	74.3
K1357-16	MGAGB3	16	1.16	9.87	7.75	75.7
K1357-17	MGAGB4	17	1.13	9.91	8.65	85.6
K1357-18	MGAGB5	18	1.14	9.92	7.68	74.5
K1357-19	MGAGB6	19	1.17	9.84	7.52	73.2
K1357-20	MGAGB7	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$