



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:25
In Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 02/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100679

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1352-01	MGAG62	1	1.14	9.96	8.28	81.0
K1352-02	MGAG63	2	1.18	9.77	7.89	78.1
K1352-03	MGAG64	3	1.19	9.96	7.75	74.8
K1352-04	MGAG65	4	1.15	9.71	7.00	68.3
K1352-05	MGAG66	5	1.17	9.94	7.45	71.6
K1352-06	MGAG67	6	1.18	9.75	7.33	71.8
K1352-07	MGAG68	7	1.14	9.5	7.03	70.5
K1352-08	MGAG69	8	1.12	9.7	7.25	71.4
K1352-09	MGAG70	9	1.18	9.98	9.05	89.4
K1352-10	MGAG70D	10	1.18	9.98	9.05	89.4
K1352-11	MGAG70S	11	1.18	9.98	9.05	89.4
K1352-12	MGAG71	12	1.16	9.94	7.83	76.0
K1352-13	MGAG72	13	1.16	9.97	7.77	75.0
K1352-14	MGAG73	14	1.13	9.62	7.93	80.1
K1352-15	MGAG74	15	1.16	9.8	7.99	79.1
K1352-16	MGAG75	16	1.17	9.95	7.79	75.4
K1352-17	MGAG76	17	1.14	9.57	7.51	75.6
K1352-18	MGAG77	18	1.18	9.94	7.76	75.1
K1352-19	MGAG78	19	1.13	9.91	8.61	85.2
K1352-20	MGAG79	20	1.17	9.64	7.64	76.4
K1352-21	MGAG80	21	1.16	9.69	7.97	79.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$