



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:35
In Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:34
Out Date: 02/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100680

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1354-01	MGAG81	1	1.1	9.54	8.31	85.4
K1354-02	MGAG82	2	1.17	9.82	8.33	82.8
K1354-03	MGAG83	3	1.18	9.62	7.9	79.6
K1354-04	MGAG84	4	1.17	9.66	8.37	84.8
K1354-05	MGAG85	5	1.14	9.77	8.44	84.6
K1354-06	MGAG86	6	1.16	9.71	8.52	86.1
K1354-07	MGAG87	7	1.11	9.94	7.97	77.7
K1354-08	MGAG88	8	1.19	9.92	8.01	78.1
K1354-09	MGAG89	9	1.18	9.97	9.08	89.9
K1354-10	MGAG90	10	1.17	9.94	8.3	81.3
K1354-11	MGAG90D	11	1.17	9.94	8.3	81.3
K1354-12	MGAG90S	12	1.17	9.94	8.3	81.3
K1354-13	MGAG91	13	1.13	9.9	8.04	78.8
K1354-14	MGAG92	14	1.14	9.98	8.12	79.0
K1354-15	MGAG93	15	1.15	9.96	7.78	75.3
K1354-16	MGAG94	16	1.19	9.97	8.99	88.8
K1354-17	MGAG95	17	1.11	9.94	9.13	90.8
K1354-18	MGAG96	18	1.18	9.93	8.75	86.5
K1354-19	MGAG97	19	1.19	9.98	7.18	68.1
K1354-20	MGAG98	20	1.17	9.66	7.76	77.6
K1354-21	MGAG99	21	1.16	9.91	8.16	80.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$