



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 02/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100699

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1383-01	A0C5B-SS1-0.0-0.5	1	1.13	9.6	8.39	85.7
K1383-02	A0C5B-SS1-4.5-5.0	2	1.15	9.7	8.63	87.5
K1383-03	A0C5B-SS2-0.0-0.5	3	1.16	9.92	8.41	82.8
K1385-01	98-SB-01-2-4	4	1.13	9.78	9.11	92.3
K1385-02	K1385-01MS	5	1.13	9.78	9.11	92.3
K1385-03	K1385-01MSD	6	1.13	9.78	9.11	92.3
K1385-04	98-SB-101-2-4	7	1.15	9.88	9.18	92.0
K1385-05	98-SB-01-5-7	8	1.12	9.97	9.32	92.7
K1385-06	98-SB-01-4-57	9	1.16	9.86	8.82	88.0
K1385-07	98-SB-01-55-57	10	1.15	9.98	7.91	76.6
K1389-01	M-4-1	11	1.13	9.57	8.8	90.9
K1389-02	M-4-2	12	1.16	9.96	9.32	92.7
K1389-03	M-4	13	1.1	9.95	9.35	93.2
K1389-05	M-6-1	14	1.15	9.98	9.28	92.1
K1389-06	M-6-2	15	1.16	9.97	9.01	89.1
K1389-07	M-6	16	1.12	9.77	8.82	89.0
K1389-09	M-5-S1	17	1.14	9.79	8.86	89.2
K1389-10	M-5-S2	18	1.15	9.75	8.81	89.1
K1389-11	M-5S	19	1.12	9.51	8.61	89.3
K1389-13	M-5-D1	20	1.13	9.88	9.00	89.9
K1389-14	M-5-D2	21	1.16	9.96	9.15	90.8
K1389-15	M-5-D	22	1.11	9.75	8.89	90.0
K1389-17	M-7-S1	23	1.15	9.71	8.6	87.0
K1389-18	M-7-S2	24	1.15	9.95	9.03	89.5
K1389-19	M-7-S	25	1.16	9.94	8.93	88.5
K1393-01	OILY-COMP-1	26	1.00	2.00	2.00	100.0
K1393-02	OILY-DEBRIS	27	1.00	2.00	2.00	100.0
K1393-03	OILY-COMP-2	28	1.00	2.00	2.00	100.0
K1394-01	ISW	29	1.00	2.00	2.00	100.0
K1394-03	ISF	30	1.00	2.00	2.00	100.0
K1396-01	23935	31	1.14	9.98	8.13	79.1
K1401-01	RBR-200050	32	1.12	9.69	8.43	85.3
K1401-03	RBR-200028	33	1.15	9.97	8.17	79.6
K1401-05	RT-3425	34	1.15	9.96	8.3	81.2
K1401-07	RT-3451	35	1.15	9.98	9.22	91.4
K1401-09	RBR-200052	36	1.16	9.57	7.8	79.0
K1403-01	S-1	37	1.15	9.55	8.77	90.7
K1403-02	S-1A	38	1.16	9.68	8.98	91.8



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 02/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 02/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100699

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$