



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 2/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 17:35

In Date: 02/08/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 02/09/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB100728

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1348-01	OK-01-020719	1	1.13	9.91	9.03	90.0
K1349-01	OR-03-020719-A	2	1.15	9.66	8.74	89.2
K1349-03	OR-03-020719-B	3	1.18	9.74	8.7	87.9
K1356-05	SU-03-020819-A	81	1.19	9.88	9.08	90.8
K1356-07	SU-03-020819-B	82	1.12	9.55	8.11	82.9
K1385-08	98-WC-S-01-2-57	4	1.12	9.71	9.35	95.8
K1390-01	SB-01 (12-13)	5	1.14	9.94	8.27	81.0
K1390-02	SB-02 (15-16)	6	1.18	9.96	8.31	81.2
K1390-03	SB-03 (12-13)	7	1.16	9.84	8.02	79.0
K1390-04	SB-04 (0-2)	8	1.18	9.95	8.38	82.1
K1390-05	SB-04 (13-14)	9	1.11	9.82	8.08	80.0
K1390-06	SB-05 (0-2)	10	1.12	9.95	7.73	74.9
K1390-07	SB-05 (8-9)	11	1.14	9.69	8.03	80.6
K1390-08	SB-06 (0-2)	12	1.17	9.87	7.36	71.1
K1390-09	SB-06 (13-14)	13	1.14	9.96	8.34	81.6
K1390-12	SB-07 (0-2)	14	1.14	9.68	7.6	75.6
K1390-13	SB-07 (7-8)	15	1.19	9.97	8.43	82.5
K1390-14	DUPE-X	16	1.14	9.98	8.01	77.7
K1390-15	DUPE-Y	17	1.16	9.94	8.74	86.3
K1404-02	AOC12-SS1-7.5-8.0	18	1.14	9.53	8.14	83.4
K1404-03	K1404-02MS	19	1.14	9.53	8.14	83.4
K1404-04	K1404-02MSD	20	1.14	9.53	8.14	83.4
K1404-05	DUP-02-06-19	21	1.11	9.82	8.26	82.1
K1404-06	AOC12-SS2-3.5-4.0	22	1.18	9.96	8.85	87.4
K1404-07	AOC12-SS3-9.0-9.5	23	1.17	9.62	8.28	84.1
K1404-08	AOC12-SS4-10.0-10.5	24	1.13	9.98	8.62	84.6
K1409-01	SB-08 (0-2)	25	1.15	9.72	8.75	88.7
K1409-02	SB-08 (14-17)	26	1.17	9.93	8.9	88.2
K1409-03	SB-09 (0-2)	27	1.1	9.57	8.67	89.4
K1409-04	SB-09 (13-16)	28	1.15	9.95	8.37	82.0
K1411-01	EIGI-RO-294	29	1.16	9.63	8.88	91.1
K1412-01	TP-18	30	1.16	9.92	8.83	87.6
K1412-03	TP-17	31	1.17	9.57	8.00	81.3
K1412-05	TP-15	32	1.19	9.79	8.54	85.5
K1413-01	DRILL-HALL-EAST	33	1.00	2.00	2.00	100.0
K1413-02	BLANK	34	1.00	2.00	2.00	100.0
K1422-01	TEST-PITS-1	35	1.15	9.93	8.64	85.3
K1422-02	TEST-PITS-1-A	36	1.15	9.93	8.64	85.3
K1422-03	TEST-PITS-1-B	37	1.15	9.93	8.64	85.3



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 2/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 17:35

In Date: 02/08/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 02/09/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB100728

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1422-04	TEST-PITS-1	61	1.15	9.93	8.64	85.3
K1422-05	CONCRETE	38	1.00	2.00	2.00	100.0
K1422-06	PAINT-CHIPS	62	1.00	2.00	2.00	100.0
K1423-01	HB-VJ210	39	1.15	9.96	8.87	87.6
K1423-02	HB-VJ210	40	1.15	9.96	8.87	87.6
K1424-01	SS2-RBR200035	41	1.13	9.93	5.84	53.5
K1424-02	SS2-RBR200035	42	1.13	9.93	5.84	53.5
K1425-01	PM-1	43	1.00	2.00	2.00	100.0
K1425-02	PM-2	44	1.00	2.00	2.00	100.0
K1425-03	PM-3	45	1.00	2.00	2.00	100.0
K1425-04	PM-4	46	1.00	2.00	2.00	100.0
K1425-05	PM-5	47	1.00	2.00	2.00	100.0
K1425-06	PM-6	48	1.00	2.00	2.00	100.0
K1425-07	PM-7	49	1.00	2.00	2.00	100.0
K1425-08	PM-8	50	1.00	2.00	2.00	100.0
K1425-09	PM-9	51	1.00	2.00	2.00	100.0
K1425-10	PM-10	52	1.00	2.00	2.00	100.0
K1425-11	PM-11	53	1.00	2.00	2.00	100.0
K1425-12	PM-12	54	1.00	2.00	2.00	100.0
K1425-13	PM-13	55	1.00	2.00	2.00	100.0
K1425-14	PM-14	56	1.00	2.00	2.00	100.0
K1425-15	PM-15	57	1.00	2.00	2.00	100.0
K1425-16	PM-16	58	1.00	2.00	2.00	100.0
K1425-17	PM-17	59	1.00	2.00	2.00	100.0
K1425-18	PM-18	60	1.00	2.00	2.00	100.0
K1425-19	PM-19	63	1.00	2.00	2.00	100.0
K1425-20	PM-20	64	1.00	2.00	2.00	100.0
K1442-01	VOC-10	65	1.13	6.75	6.22	90.6
K1442-02	VOC-11	66	1.13	5.88	5.23	86.3
K1442-03	VOC-12	67	1.15	6.47	5.74	86.3
K1442-04	COMP-10	68	1.19	9.74	8.66	87.4
K1442-05	COMP-11	69	1.19	9.96	8.84	87.2
K1442-06	COMP-12	70	1.17	9.91	8.72	86.4
K1445-05	SVOC-GPC-BLANK	71	1.00	2.00	2.00	100.0
K1445-06	PEST-GPC-BLANK	72	1.00	2.00	2.00	100.0
K1445-07	PEST-GPC-BLANK-SPIKE	73	1.00	2.00	2.00	100.0
K1445-08	PCB-GPC-BLANK	74	1.00	2.00	2.00	100.0
K1445-09	PCB-GPC-BLANK-SPIKE	75	1.00	2.00	2.00	100.0
K1445-10	SVOC-GPC2-BLANK	76	1.00	2.00	2.00	100.0



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 2/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 17:35

In Date: 02/08/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 02/09/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB100728

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1445-11	PEST-GPC2-BLANK	77	1.00	2.00	2.00	100.0
K1445-12	PEST-GPC2-BLANK-SPIKE	78	1.00	2.00	2.00	100.0
K1445-13	PCB-GPC2-BLANK	79	1.00	2.00	2.00	100.0
K1445-14	PCB-GPC2-BLANK-SPIKE	80	1.00	2.00	2.00	100.0
K1450-01	SHORE-RD-BOX-CULVERT-TP	83	1.14	9.76	8.07	80.4
K1450-03	HOSPITAL-BOX-CULVERT-TP	84	1.14	9.95	8.78	86.7
K1452-01	TEST-PIT-2	85	1.13	9.84	8.79	87.9
K1452-02	TEST-PIT-2-A	86	1.13	9.84	8.79	87.9
K1452-03	TEST-PIT-2-B	87	1.13	9.84	8.79	87.9
K1452-04	TEST-PIT-2	88	1.13	9.84	8.79	87.9
K1453-01	M-8-S1	89	1.14	9.76	8.85	89.4
K1453-02	M-8-S2	90	1.16	9.72	8.67	87.7
K1453-03	M-8-S	91	1.15	9.84	8.86	88.7
K1453-05	M-8-D1	92	1.13	9.6	8.5	87.0
K1453-06	M-8-D2	93	1.1	9.96	8.96	88.7
K1453-07	M-8-D	94	1.18	9.98	9.12	90.2
K1453-09	M-9-S1	95	1.15	9.61	8.65	88.7
K1453-10	M-9-S2	96	1.17	9.66	8.57	87.2
K1453-11	M-9-S	97	1.17	9.92	9.00	89.5
K1453-13	M-9-D1	98	1.16	9.97	8.92	88.1
K1453-14	M-9-D2	99	1.12	9.83	8.02	79.2
K1453-15	M-9-D	100	1.13	9.85	8.00	78.8
K1454-01	TP-1	101	1.14	9.96	8.92	88.2
K1454-02	TP-1-A	102	1.14	9.96	8.92	88.2
K1454-03	TP-1-B	103	1.14	9.96	8.92	88.2
K1454-04	TP-1	104	1.14	9.96	8.92	88.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$