



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 2/12/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:10
In Date: 02/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:10
Out Date: 02/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100761

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1430-01	MH0YC7	1	1.00	2.00	2.00	100.0
K1430-02	MH0YC8	2	1.00	2.00	2.00	100.0
K1430-03	MH0YC9	3	1.00	2.00	2.00	100.0
K1430-04	MH0YD0	4	1.00	2.00	2.00	100.0
K1430-05	MH0YD1	5	1.00	2.00	2.00	100.0
K1430-06	MH0YD2	6	1.00	2.00	2.00	100.0
K1430-07	MH0YD3	7	1.00	2.00	2.00	100.0
K1430-08	MH0YD4	8	1.00	2.00	2.00	100.0
K1430-09	MH0YD5	9	1.00	2.00	2.00	100.0
K1430-10	MH0YD6	10	1.00	2.00	2.00	100.0
K1430-11	MH0YD6D	11	1.00	2.00	2.00	100.0
K1430-12	MH0YD6S	12	1.00	2.00	2.00	100.0
K1430-13	MH0YD7	13	1.00	2.00	2.00	100.0
K1430-14	MH0YD8	14	1.00	2.00	2.00	100.0
K1430-15	MH0YJ3	15	1.00	2.00	2.00	100.0
K1430-16	MH0YJ4	16	1.00	2.00	2.00	100.0
K1430-17	MH0YJ5	17	1.00	2.00	2.00	100.0
K1430-18	MH0YJ6	18	1.00	2.00	2.00	100.0
K1430-19	MH0YJ7	19	1.00	2.00	2.00	100.0
K1430-20	MH0YJ8	20	1.00	2.00	2.00	100.0
K1430-21	MH0YJ9	21	1.00	2.00	2.00	100.0
K1430-22	MH0YK0	22	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$