



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 2/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 02/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 02/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100841

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1346-07	CL-01-021419-A	1	1.13	9.78	8.9	89.8
K1346-09	CL-01-021419-B	2	1.15	9.6	8.52	87.2
K1346-11	CL-01-021419-C	3	1.14	9.96	9.19	91.3
K1346-13	CL-01-021419-D	4	1.17	9.77	8.98	90.8
K1346-15	CL-01-021419-E	5	1.19	9.98	9.26	91.8
K1358-03	TR-04-021419	6	1.15	9.67	8.84	90.3
K1484-01	M-2-S1	7	1.14	9.88	9.1	91.1
K1484-02	M-2-S2	8	1.17	9.94	9.23	91.9
K1484-03	M-2-S	9	1.16	9.93	9.21	91.8
K1484-05	M-2-D1	10	1.13	9.7	8.7	88.3
K1484-06	M-2-D2	11	1.14	9.94	8.76	86.6
K1484-07	M-2-D	12	1.16	9.91	8.89	88.3
K1484-09	M-3-S1	13	1.13	9.97	9.31	92.5
K1484-10	M-3-S2	14	1.18	9.63	8.39	85.3
K1484-11	M-3-S	15	1.19	9.94	8.78	86.7
K1484-13	M-3-D1	16	1.17	9.97	9.11	90.2
K1484-14	M-3-D2	17	1.17	9.69	8.78	89.3
K1484-15	M-3-D	18	1.16	9.67	8.78	89.5
K1486-01	VNJ-227	19	1.18	9.98	8.53	83.5
K1486-02	VNJ-227	20	1.18	9.98	8.53	83.5
K1487-01	RHR-1	21	1.00	2.00	2.00	100.0
K1490-01	TP-1-214	22	1.16	9.71	8.28	83.3
K1490-02	TP-1-214-A	23	1.16	9.71	8.28	83.3
K1490-03	TP-1-214-B	24	1.16	9.71	8.28	83.3
K1490-04	TP-1-214	25	1.16	9.71	8.28	83.3
K1498-01	A0C14-02-B3 (6-12)	26	1.13	9.55	8.8	91.1
K1498-02	A0C14-02-E (6-12)	27	1.17	9.75	8.1	80.8
K1498-03	A0C14-02-B4 (0-6)	28	1.19	9.54	8.66	89.5
K1498-04	A0C14-02-E1 (0-6)	29	1.19	9.77	9.02	91.3
K1498-05	A0C14-02-C4 (0-6)	30	1.11	9.57	9.09	94.3
K1498-06	A0C14-02-F (0-6)	31	1.16	9.68	8.27	83.5
K1498-07	A0C14-02-D2 (0-6)	32	1.15	9.95	7.93	77.0
K1498-08	A0C2-02-F1 (0-6)	33	1.14	9.94	8.75	86.5
K1498-09	A0C2-04-F1 (6-12)	34	1.13	9.93	8.59	84.8
K1498-10	A0C2-04-E1 (0-6)	35	1.11	9.91	7.88	76.9
K1498-11	A0C2-04-E1 (6-12)	36	1.19	9.97	8.77	86.3
K1498-12	A0C2-04-C1 (12-18)	37	1.12	9.93	8.69	85.9
K1498-13	A0C2-04-G1 (0-6)	38	1.19	9.95	9.11	90.4
K1498-14	A0C2-04-G1 (6-12)	39	1.14	9.98	9.37	93.1



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 2/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 02/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 02/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100841

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1500-01	G-021419	40	1.14	9.66	7.95	79.9
K1500-02	G-021419	41	1.14	9.66	7.95	79.9
K1500-05	G-021419	42	1.14	9.66	7.95	79.9
K1501-01	CH-021419	43	1.17	9.96	8.61	84.6
K1501-02	CH-021419	44	1.17	9.96	8.61	84.6
K1501-05	CH-021419	45	1.17	9.96	8.61	84.6
K1502-01	RO-1-4	46	1.19	9.88	7.99	78.3
K1502-02	RO-1-4	47	1.19	9.88	7.99	78.3
K1502-04	DRUM-A-1-5	48	1.00	2.00	2.00	100.0
K1502-05	DRUM-B-1-5	49	1.00	2.00	2.00	100.0
K1502-06	DRUM-C-1-5	50	1.00	2.00	2.00	100.0
K1502-07	DRUM-D-1-5	51	1.00	2.00	2.00	100.0
K1502-08	DRUM-E-1-5	52	1.00	2.00	2.00	100.0
K1502-09	DRUM-F-1-5	53	1.00	2.00	2.00	100.0
K1502-10	DRUM-G-1-5	54	1.00	2.00	2.00	100.0
K1502-11	DRUM-H-1-5	55	1.00	2.00	2.00	100.0
K1502-12	DRUM-I-1-5	56	1.00	2.00	2.00	100.0
K1502-13	DRUM-J-1-5	57	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$