



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 2/23/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:30
In Date: 02/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:32
Out Date: 02/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100962

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1573-01	MH0YM9	1	1.00	2.00	2.00	100.0
K1573-02	MH0YN0	2	1.00	2.00	2.00	100.0
K1573-03	MH0YN1	3	1.00	2.00	2.00	100.0
K1573-04	MH0YN2	4	1.00	2.00	2.00	100.0
K1573-05	MH0YN2D	5	1.00	2.00	2.00	100.0
K1573-06	MH0YN2S	6	1.00	2.00	2.00	100.0
K1573-07	MH0YN3	7	1.00	2.00	2.00	100.0
K1573-08	MH0YN4	8	1.00	2.00	2.00	100.0
K1573-09	MH0YN5	9	1.00	2.00	2.00	100.0
K1573-10	MH0YN6	10	1.00	2.00	2.00	100.0
K1573-11	MH0YN7	11	1.00	2.00	2.00	100.0
K1573-12	MH0YN8	12	1.00	2.00	2.00	100.0
K1573-13	MH0YN9	13	1.00	2.00	2.00	100.0
K1573-14	MH0YP0	14	1.00	2.00	2.00	100.0
K1573-15	MH0YP1	15	1.00	2.00	2.00	100.0
K1573-16	MH0YP2	16	1.00	2.00	2.00	100.0
K1573-17	MH0YP3	17	1.00	2.00	2.00	100.0
K1573-18	MH0YP4	18	1.00	2.00	2.00	100.0
K1573-19	MH0YP5	19	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$