



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:00
In Date: 02/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 02/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100963

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1574-01	MH0YP6	1	1.00	2.00	2.00	100.0
K1574-02	MH0YP7	2	1.00	2.00	2.00	100.0
K1574-03	MH0YP8	3	1.00	2.00	2.00	100.0
K1574-04	MH0YP9	4	1.00	2.00	2.00	100.0
K1574-05	MH0YP9D	5	1.00	2.00	2.00	100.0
K1574-06	MH0YP9S	6	1.00	2.00	2.00	100.0
K1574-07	MH0YQ0	7	1.00	2.00	2.00	100.0
K1574-08	MH0YQ1	8	1.00	2.00	2.00	100.0
K1574-09	MH0YQ2	9	1.00	2.00	2.00	100.0
K1574-10	MH0YQ3	10	1.00	2.00	2.00	100.0
K1574-11	MH0YQ4	11	1.00	2.00	2.00	100.0
K1574-12	MH0YQ5	12	1.00	2.00	2.00	100.0
K1574-13	MH0YQ6	13	1.00	2.00	2.00	100.0
K1574-14	MH0YQ7	14	1.00	2.00	2.00	100.0
K1574-15	MH0YQ8	15	1.00	2.00	2.00	100.0
K1574-16	MH0YQ9	16	1.00	2.00	2.00	100.0
K1574-17	MH0YR0	17	1.00	2.00	2.00	100.0
K1574-18	MH0YR1	18	1.00	2.00	2.00	100.0
K1574-19	MH0YR2	19	1.00	2.00	2.00	100.0
K1574-20	MH0YR3	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$