



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 2/23/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:40
In Date: 02/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 02/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB100973

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1579-01	DB5Z1	1	1.00	2.00	2.00	100.0
K1579-02	DB5Z8	2	1.14	9.78	7.51	73.7
K1579-03	DB600	3	1.15	9.96	4.04	32.8
K1579-04	DB601	4	1.13	9.55	5.1	47.1
K1579-05	DB602	5	1.12	9.96	8.03	78.2
K1579-06	DB603	6	1.15	9.91	8.11	79.5
K1579-07	DB604	7	1.15	9.73	7.08	69.1
K1579-08	DB605	8	1.14	9.97	9.04	89.5
K1579-09	DB606	9	1.17	9.93	9.00	89.4
K1579-10	DB607	10	1.17	9.96	8.56	84.1
K1579-11	DB608	11	1.16	9.91	8.41	82.9
K1579-12	DB609	12	1.19	9.92	7.53	72.6
K1579-13	DB609MS	13	1.19	9.92	7.53	72.6
K1579-14	DB609MSD	14	1.19	9.92	7.53	72.6
K1579-15	DB610	15	1.17	9.58	8.33	85.1
K1579-16	DB611	16	1.17	9.64	8.69	88.8
K1579-17	DB612	17	1.15	9.8	9.72	99.1
K1579-18	VHBLK01	18	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$