



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/1/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 02/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 03/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101090

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1655-01	A0C12-1-1.0-1.5	1	1.14	9.96	9.00	89.1
K1655-02	DUP-2	2	1.18	9.94	9.02	89.5
K1655-03	A0C1-1-17.5-18.0	3	1.15	9.78	8.76	88.2
K1655-10	A0C12-2-1.5-2.0	4	1.17	9.6	8.22	83.6
K1655-11	K1655-10MS	5	1.17	9.6	8.22	83.6
K1655-12	K1655-10MSD	6	1.17	9.6	8.22	83.6
K1655-13	A0C1-2-5.5-6.0	7	1.17	9.79	9.13	92.3
K1655-14	A0C1-2-14.5-15.0	8	1.14	9.97	8.92	88.1
K1665-01	200051	9	1.14	9.81	7.5	73.4
K1665-03	200029	10	1.14	9.95	8.62	84.9
K1665-05	200053	11	1.18	9.86	7.84	76.7
K1665-07	200031	12	1.13	9.94	7.86	76.4
K1666-01	PM-A	13	1.00	2.00	2.00	100.0
K1666-02	PM-B	14	1.00	2.00	2.00	100.0
K1666-03	PM-C	15	1.00	2.00	2.00	100.0
K1666-04	PM-D	16	1.00	2.00	2.00	100.0
K1666-05	PM-E	17	1.00	2.00	2.00	100.0
K1666-06	PM-F	18	1.00	2.00	2.00	100.0
K1666-07	PM-G	19	1.00	2.00	2.00	100.0
K1666-08	PM-H	20	1.00	2.00	2.00	100.0
K1669-01	FK-GF-02-031-A	21	1.15	9.94	8.66	85.4
K1669-02	FK-GF-02-031-A	22	1.17	9.54	8.28	84.9
K1669-03	FK-GF-02-031-B	23	1.19	9.61	8.43	86.0
K1669-04	FK-GF-02-031-B	24	1.14	9.94	8.83	87.4
K1670-01	FK-GF-02-032-A	25	1.14	9.96	8.7	85.7
K1670-02	FK-GF-02-032-A	26	1.14	9.98	8.78	86.4
K1670-03	FK-GF-02-032-B	27	1.12	9.64	8.37	85.1
K1670-04	FK-GF-02-032-B	28	1.15	9.65	8.29	84.0
K1671-01	FK-BPM-04-002-ABCD	29	1.19	9.63	8.62	88.0
K1671-02	FK-BPM-04-002-ABCD	30	1.14	9.97	9.28	92.2
K1672-01	FK-GF-02-033-A	31	1.18	9.97	8.97	88.6
K1672-02	FK-GF-02-033-A	32	1.13	9.98	8.66	85.1
K1672-03	FK-GF-02-033-B	33	1.15	9.69	8.31	83.8
K1672-04	FK-GF-02-033-B	34	1.1	9.97	8.6	84.6
K1675-01	Z-3-11-A	35	1.15	9.74	9.25	94.3
K1675-02	Z-3-11-A-1	36	1.17	9.92	9.44	94.5
K1675-03	Z-3-11-A-2	37	1.17	9.92	9.44	94.5
K1675-04	Z-3-11-A	38	1.15	9.74	9.25	94.3
K1675-05	Z-3-11-B	39	1.16	9.93	9.27	92.5



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/1/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 02/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 03/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101090

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1675-06	Z-3-11-B-1	40	1.19	9.74	9.31	95.0
K1675-07	Z-3-11-B-2	41	1.19	9.74	9.31	95.0
K1675-08	Z-3-11-B	42	1.16	9.93	9.27	92.5
K1675-09	Z-3-1-A	43	1.1	9.68	8.81	89.9
K1675-10	Z-3-1-A	44	1.1	9.68	8.81	89.9
K1675-11	Z-3-1-B	45	1.17	9.96	9.1	90.2
K1675-12	Z-3-1-B	46	1.17	9.96	9.1	90.2
K1675-13	Z-3-2-A	47	1.16	9.94	9.26	92.3
K1675-15	Z-3-2-B	48	1.16	9.96	9.44	94.1
K1675-16	Z-3-2-B	49	1.16	9.96	9.44	94.1
K1675-17	Z-4-A-B	50	1.18	9.98	9.11	90.1
K1675-18	Z-4-A-B	51	1.18	9.98	9.11	90.1
K1675-19	Z-4-1-A-1	52	1.13	9.67	8.85	90.4
K1675-20	Z-4-1-A-2	53	1.13	9.67	8.85	90.4
K1675-21	Z-4-1-B-1	54	1.19	9.81	8.82	88.5
K1675-22	Z-4-1-B-2	55	1.19	9.81	8.82	88.5
K1675-23	Z-3-12-A	56	1.14	9.96	9.13	90.6
K1675-24	Z-3-12-A-1	57	1.16	9.85	9.15	91.9
K1675-25	Z-3-12-A-2	58	1.16	9.85	9.15	91.9
K1675-26	Z-3-12-A	59	1.14	9.96	9.13	90.6
K1675-27	Z-3-12-B	60	1.14	9.76	8.93	90.4
K1675-28	Z-3-12-B	61	1.14	9.76	8.93	90.4
K1675-29	Z-3-13-A	62	1.15	9.95	9.22	91.7
K1675-30	Z-3-13-A-1	63	1.17	9.98	9.31	92.4
K1675-31	Z-3-13-A-2	64	1.17	9.98	9.31	92.4
K1675-32	Z-3-13-A	65	1.15	9.95	9.22	91.7
K1675-33	Z-3-13-B	66	1.17	9.63	8.52	86.9
K1675-34	Z-3-13-B-1	67	1.2	9.72	8.87	90.0
K1675-35	Z-3-13-B-2	68	1.2	9.72	8.87	90.0
K1675-36	Z-3-13-B	69	1.17	9.63	8.52	86.9
K1675-37	Z-3-14-A	70	1.16	9.82	9.09	91.6
K1675-38	Z-3-14-A-1	71	1.11	9.6	9.08	93.9
K1675-39	Z-3-14-A-2	72	1.11	9.6	9.08	93.9
K1675-40	Z-3-14-A	73	1.16	9.82	9.09	91.6
K1675-41	Z-3-14-B	74	1.17	9.75	9.12	92.7
K1675-42	Z-3-14-B-1	75	1.13	9.72	8.74	88.6
K1675-43	Z-3-14-B-2	76	1.13	9.72	8.74	88.6
K1675-44	Z-3-14-B	77	1.17	9.75	9.12	92.7
K1675-46	Z-3-12-B	78	1.1	9.79	9.01	91.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/1/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:30
In Date: 02/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 03/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101090

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1675-47	Z-3-12-B	79	1.1	9.79	9.01	91.0
K1679-01	390207017.1	80	1.16	9.96	1.43	3.1
K1679-02	390207019	81	1.12	9.78	1.6	5.5
K1679-03	390207743	82	1.16	9.5	2.52	16.3
K1679-04	390208060	83	1.18	9.81	5.61	51.3
K1679-05	390208061	84	1.17	9.98	1.57	4.5
K1680-01	390208085	85	1.11	9.79	1.23	1.4
K1681-01	390205027	86	1.16	9.86	2.96	20.7
K1681-02	390206772	87	1.17	9.96	1.61	5.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$