



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 3/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 03/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:20
Out Date: 03/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101256

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1814-01	MH0YY8	1	1.00	2.00	2.00	100.0
K1814-02	MH0YY9	2	1.00	2.00	2.00	100.0
K1814-03	MH0YZ0	3	1.00	2.00	2.00	100.0
K1814-04	MH0YZ1	4	1.00	2.00	2.00	100.0
K1814-05	MH0YZ2	5	1.00	2.00	2.00	100.0
K1814-06	MH0YZ3	6	1.00	2.00	2.00	100.0
K1814-07	MH0YZ4	7	1.00	2.00	2.00	100.0
K1814-08	MH0YZ5	8	1.00	2.00	2.00	100.0
K1814-09	MH0YZ6	9	1.00	2.00	2.00	100.0
K1814-10	MH0YZ7	10	1.00	2.00	2.00	100.0
K1814-11	MH0YZ8	11	1.00	2.00	2.00	100.0
K1814-12	MH0YZ9	12	1.00	2.00	2.00	100.0
K1814-13	MH0Z00	13	1.00	2.00	2.00	100.0
K1814-14	MH0Z01	14	1.00	2.00	2.00	100.0
K1814-15	MH0Z02	15	1.00	2.00	2.00	100.0
K1814-16	MH0Z03	16	1.00	2.00	2.00	100.0
K1814-17	MH0Z04	17	1.00	2.00	2.00	100.0
K1814-18	MH0Z04D	18	1.00	2.00	2.00	100.0
K1814-19	MH0Z04S	19	1.00	2.00	2.00	100.0
K1814-20	MH0Z05	20	1.00	2.00	2.00	100.0
K1814-21	MH0Z06	21	1.00	2.00	2.00	100.0
K1814-22	MH0Z07	22	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$